

## Product datasheet for **RR212231L3V**

### Ramp2 (NM\_031646) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ramp2 (NM_031646) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ramp2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031646
ORF Size:	546 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212231).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031646.1</a> , <a href="#">NP_113834.1</a>
RefSeq Size:	751 bp
RefSeq ORF:	549 bp
Locus ID:	58966
UniProt ID:	<a href="#">Q9JHJ1</a>
Cytogenetics:	10q31
Gene Summary:	transports calcitonin receptor-like receptor (Calcr) to the cell membrane and determine its glycosylation state and ligand specificity as a functional adrenomedullin receptor [RGD, Feb 2006]



[View online »](#)