

Product datasheet for **RR212236L4V**

Ptk2 (NM_013081) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ptk2 (NM_013081) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ptk2
Synonyms:	FAK; FRNK; p125FAK
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013081
ORF Size:	3165 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212236).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013081.2 , NP_037213.1
RefSeq Size:	4559 bp
RefSeq ORF:	3168 bp
Locus ID:	25614
UniProt ID:	O35346
Cytogenetics:	7q34
Gene Summary:	tyrosine kinase that is enriched in focal adhesions; may regulate mechanical signal transduction in cardiac myocytes [RGD, Feb 2006]



[View online »](#)