

## Product datasheet for **RR212241L3V**

### Sycp3 (NM\_013041) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sycp3 (NM_013041) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sycp3
Synonyms:	RNASCP3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013041
ORF Size:	771 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212241).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013041.1</a> , <a href="#">NP_037173.1</a>
RefSeq Size:	1132 bp
RefSeq ORF:	774 bp
Locus ID:	25561
UniProt ID:	<a href="#">Q63520</a>
Cytogenetics:	7q13
Gene Summary:	component of synptonemal complexes (SCs), formed between homologous chromosomes during meiotic prophase; may play a role in spermatogenesis [RGD, Feb 2006]



[View online »](#)