

Product datasheet for **RR212248L4V**

Lipf (NM_017341) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lipf (NM_017341) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lipf
Synonyms:	RNLIP; Rnlp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017341
ORF Size:	1185 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212248).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017341.2 , NP_059037.1
RefSeq Size:	1377 bp
RefSeq ORF:	1188 bp
Locus ID:	50682
UniProt ID:	P04634
Cytogenetics:	1q52
Gene Summary:	gastric enzyme; metabolizes lipids and triacylglycerols; may be auxiliary enzyme for fat digestion and absorption [RGD, Feb 2006]



[View online »](#)