

## Product datasheet for **RR212251L3V**

### Gja3 (NM\_024376) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gja3 (NM_024376) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gja3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024376
ORF Size:	1248 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212251).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024376.2</a> , <a href="#">NP_077352.1</a>
RefSeq Size:	1590 bp
RefSeq ORF:	1251 bp
Locus ID:	79217
UniProt ID:	<a href="#">P29414</a>
Cytogenetics:	15p12
Gene Summary:	lens gap junction protein; oligomerizes into intercellular channels and functions to exchange low molecular weight metabolites and ions between adjacent cells; induces voltage-gated currents in nonjunctional plasma membrane [RGD, Feb 2006]



[View online »](#)