

Product datasheet for **RR212261**

Hmga2 (NM_032070) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hmga2 (NM_032070) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Hmga2
Synonyms: Hmgic
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR212261 representing NM_032070
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCGCACGGGTGAGGGCGCCGGGCAGCCGTCCACATCAGCCAGGGACAACCTGCCACCCCGGGC
CACAGAAGCGAGGACGCGGCCGACCCAGGAAGCAGCCGCAAGAGCCAACCTGTGAGCCCTCTCCTAAGAG
ACCCAGAGGAAGACCCAAAGGCAGCAAGAACAAGAGCCCTCGAAAGCAGCTCAAAGAAAGCAGAGACC
ATTGGAGAAAAACGGCCAAGAGGCAGACCCAGGAAATGGCCACAACAAATTGTTCAAGAGCTGCTCAGG
AGACTGAAGAGACATCTCGCAAGAGTCCGAGAGGAAGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR212261 representing NM_032070
Red=Cloning site Green=Tags(s)

MSARGEGAGQPSTSAQQQPATPAPQKRGRGRPRKQPQEPTCEPSPKRPRGRPKGSKNKSPSKAAQKKAET
IGEKPRGRPRKWPQQIVQKPAQETEETSSQESAEED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites: Sgfl-MluI



[View online »](#)


Cloning Scheme:

Cloning sites used for ORF Shutting:



* The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_032070

ORF Size: 321 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_032070.1 , NP_114459.1
RefSeq Size:	802 bp
RefSeq ORF:	324 bp
Locus ID:	84017
Cytogenetics:	7q22
MW:	11.7 kDa
Gene Summary:	a non-histone chromosomal architecture protein [RGD, Feb 2006]