

Product datasheet for **RR212267L4V**

Gtf2ird1 (NM_001001504) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gtf2ird1 (NM_001001504) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gtf2ird1
Synonyms:	Gtf3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001001504
ORF Size:	2829 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212267).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001001504.2 , NP_001001504.1
RefSeq Size:	3110 bp
RefSeq ORF:	2832 bp
Locus ID:	246770
Cytogenetics:	12q12
Gene Summary:	may act as a muscle specific transcription factor; human homolog may be associated with the Williams-Beuren microdeletion syndrome [RGD, Feb 2006]



[View online »](#)