

## Product datasheet for **RR212278L3V**

### **MyI9 (NM\_001100885) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	MyI9 (NM_001100885) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	MyI9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001100885
ORF Size:	513 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212278).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001100885.1</a> , <a href="#">NP_001094355.1</a>
RefSeq Size:	1008 bp
RefSeq ORF:	516 bp
Locus ID:	296313
UniProt ID:	<a href="#">Q64122</a>
Cytogenetics:	3q42
Gene Summary:	Myosin regulatory subunit that plays an important role in regulation of both smooth muscle and nonmuscle cell contractile activity via its phosphorylation. Implicated in cytokinesis, receptor capping, and cell locomotion.[UniProtKB/Swiss-Prot Function]



[View online »](#)