

Product datasheet for **RR212294**

Klf3 (NM_001105742) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Klf3 (NM_001105742) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Klf3
Synonyms:	Bklf; Tef-2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR212294 representing NM_001105742 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTCATGTTTGACCCAGTCCCTGTCAAGCAAGAGGCCATGGACCCCGTCACGGTGTCTACCCATCAA
ATTACATAGAATCAATGAAGCCCAACAATATGGGGTCATTTACTCTACACCATTACCTGACAAGTTCTT
CCAGACCCAGAAGGCCCTCGCTCAGGGATACAGATGGAGCCTGTGGACCTCACTGTGAACAAGCGGGC
TCCCCGCCTTCTGCTGCTGGTTCTCCTTCTCTGAAATTTCCATCACACCGGAGAGCATCACCGGGC
TCAGCATGCCTTCTCCAGCCACCCATTAAGAAATACTACCCCCCTCTCTGTTGTGCAACCTTCGG
AGTACCACTGTCCATGCCACCTGTGATGGCCGCTGCACTGACCAGACATGGGATCCGGAGCCCGGGCATC
CTCCCTGTATCCAGCCCGTCTGGTACAACCTGTCCCTTCAATGTACACCAGCCACCTGCAGCAGCCTC
TCATGGTCTCCTTGTGAGAAGAAATGGACAATCAAACAGCGGCATGCCAGTACCCGTAATTGAATCCTA
TGAGAAGCCCATATTACAGAAAAAATTAATAAGAACCTGGAATCGAACCACAGAGGACTGACTATTAC
CCTGAAGAAATGTCACCCCTTTAATGAACCTAGTGTCCCCCGCAAGCGTTGTTGCAAGAGAACCACC
CTTCCGTATTGTGACGCCGGGAAGAGACCTCTACCTGTGGAATCCCGGATACCCAAAGGAAGCGCAG
GATACACAGATGCGACTATGATGGATGCAATAAGGTCTACACTAAGAGCTCGCACTTGAAAGCACACAGA
AGAACCATACAGGGGAAAAGCCGTACAATGCACCTGGGAAGGCTGCACGTGGAAGTTTCTCGTCTG
ATGAACTGACGAGACATTTCCGAAAACACTGGAATCAAACCTTTCCAGTGTCCAGACTGTGACCGAAG
CTTCTCTCGCTCCGACCCTTGCCCTACATAGGAAACGCCACATGCTAGTC

ACGGTACGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR212294 representing NM_001105742
 Red=Cloning site Green=Tags(s)

MLMFDVPVKQEAMDPVTVSYPSNYIESMKPNKYGVIYSTPLPDKFFQTPEGLAHGIQMEPVDLTVNKRGSPPSAAGSPSSLKFPSHRRASPLSMPSSSPPIKKYSPSPGVQPFVPLSMPPVMAAALTRHGIRSPGILPVIQPVVVQPVPFMYTSHLQQPLMVSLSSEMDNSNSGMPVPIVIESYEKPIQKKIKIEPGIEPQRDYYPEEMSPPLMNSVSPQALLQENHPSVIVQPGKRPLPVESPDTRKRRRIHRCDYDGCNKVYTKSSHLKAHRTHHTGEKPYKCTWEGCTWKFARSDDELTRHFRKHTGIKPFQCPDCDRSFSRSDHLALHRKRHMLV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001105742

ORF Size: 1032 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001105742.1 , NP_001099212.1
RefSeq Size:	4988 bp
RefSeq ORF:	1035 bp
Locus ID:	114845
Cytogenetics:	14p11
MW:	38.6 kDa
Gene Summary:	human homolog is a transcription factor; involved in haematopoietic differentiation [RGD, Feb 2006]