

Product datasheet for **RR212297**

Srpk1 (NM_001025726) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Srpk1 (NM_001025726) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Srpk1
Synonyms:	MGC116186
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:


>RR212297 representing NM_001025726
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGCGGAAAGTGCTCGCGCTCCAGGCCGAAAGAAAAGGACCAAGGCCAAGAAGGACAAAGCCCAA
GGAAGCCTGAAACCCAGCACCGAGGCTCTGCACCCCACTCCGAGAGCGACATCCCGAGCAGGAGGAGGA
GATTCTGGGGTCTGATGATGATGAACAGGAAGACCCTAACGACTACTGTAAGGAGGTTATCATCTTGTG
AAAATTGGAGATCTATTTAACGGGCGATACCATGTGATTCGAAAATTAGGCTGGGGACACTTTTCCACAG
TATGGTTATCATGGGATATTCAGGGAAAGAAGTTGTAGCAATGAAAGTAGTTAAAAGTGCCGAACATTA
CACAGAGACAGCACTGGATGAGATCCGGCTCCTGAAATCGGTTTCGCAATTCAGACCCGAATGATCCAAAT
GGAGAAATGGTCGTTCACTACTAGATGACTTTAAAATATCAGGCGTTAATGGAACACATATCTGCATGG
TGTTTGAAGTTTAGGACATCATCTACTCAAGTGGATCATCAAGTCCAATATCAGGGGCTCCGCTGCC
TTGTGTCAAAAAATTATTCAGCAAGTGTCCAGGGTCTGGATTACTTACACCAAGTGCCGGATCATC
CACACTGACATCAAGCCGAAAACATCCTGTTGTCCGTGAACGAGCAGTACATTCGGAGGCTGGCTGCAG
AGGCCACAGAGTGGCAACGGTCTGGAGCTCCCCACCATCGGGTCTGCAGTCAGCACTGCTCCACAGCC
TAAGCCAGCTGACAAAATGTCAAAGAATAAGAAGAAGAAATTGAAGAAGAAGCAGAAGCGCCAGGAGAA
TTACTAGAGAAGCGAATGCAGGAAATGAGGAAATGGAGAAAGAGTCGGGCCCTGGGCAAAAAAGACCTA
ACAAGCAAAGAAGATCAGAGAGTCTGTTGAAAGACCCTTGACAGAGAACCCACCTAATAAAATGACCCA
AGAGAACTTGAAGAGTCGAGTTCATTGGCCAGGATCAGACACTCACGGAGCATGGTGGAGAGGGTGGT
GCGCCAGAAATTAATTGCAATGGAGTGAATGGAATCGTTAATTACACTGCAAACAGTAACAATGAAACCT
TGAGGCTTAAAGAGGACCTACATAATGCCAATGACTGTGATGGCCGCACTCTGAAGCAGGAACCTAGCTT
CCTAAGCTCTCCAAATGGAGACAGTAGTCCGTCTCCAGAAACAGACTCTTGACACACCTTCGGGCTCGGAG
GTGCCCGAGACCATGGCGTGCCAGGCCTCTGTAGAGCGGTCACTCACCCACCCAACAGGACATCACCC
ACCTGCAGGAGCGCATCCGGGCACATACGCCTTCTGGGGATGAGCAGGAGCGGGAACACAATGGAGCGCT
GGACAGCAAAGGAAAATCTCTGCTGAAAATTTCTTATTAATCCCCTGGAGCAAAAAATGCAGAAAAG
CTCAAAGTGAAGATTGCTGATCTTGGGAATGCCTGCTGGGTGCACAAGCATTTACCCGAAGACATCCAGA
CCCGGCAGTACCGGCTTTGGAGGTCCTGATAGGATCTGGCTACAACACCCCTGCTGACATCTGGAGCAC
AGCATGCATGGCCTTTGAACTAGCCACAGGGGACTATCTGTTGAGCCTCATTGGGGGAGGATTACACG
CGAGATGAAGACCACATCGCCCTGATCATAGAGCTCCTGGGAAGGTGCCTCGAAAGCTCATTGTAGCAG
GGAATACTCCAAGGAATTTTACCAAAAAAGGTGACCTGAAACACATACCAAGTTGAAGCCTTGGGG
TCTTCTGGAGGTTCTGGTGGAGAAGTATGAGTGGCCTCAGGAGGAGGCTGCTGGCTTACAGATTTCTG
CTGCCATGTTGGAGCTCATGCCTGAGAAGAGAGCCACTGCCGCTGAGTGTCTCCGGCATCTTGGCTGA
ACTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:



ACCN: NM_001025726

ORF Size: 1965 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001025726.1](#), [NP_001020897.1](#)

RefSeq Size: 2374 bp
RefSeq ORF: 1968 bp
Locus ID: 361811
Cytogenetics: 20p12
MW: 73.9 kDa