

## Product datasheet for **RR212319**

### Rap1b (NM\_134346) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rap1b (NM\_134346) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Rap1b  
**Synonyms:** MGC93206  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR212319 representing NM\_134346  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCGTGAATATAAGCTAGTCGTTCTTGGCTCAGGAGGTGTTGGGAAGTCGGCTCTGACTGTACAGTTT  
TTCAAGGAATTTTGTGAAAAATATGATCCTACGATAGAAGATTCTTATAGAAAGCAAGTTGAAGTAGA  
TGCACAGCAGTGTATGCTTGAATCCTGGATACTGCAGGAACGGAGCAGTTTACAGCCATGAGAGATCTG  
TACATGAAGAACGGACAAGGCTTTGCGTTGGTTTACTCCATCACAGCACAGTCGACATTTAATGACTTAC  
AGGATCTAAGAGAGCAGATTCTTCGGGTTAAAGACACTGACGATGTCCCAATGATTCTGGTTGGTAACAA  
ATGCGACTGGAAGACGAAAGAGTTGTCCGGAAGGAACAAGGTCAGAACCTGGCAAGACAGTGGAGCAAC  
TGCGCTTTCTTAGAGTCTCCGCAAGTCCAAGATAAATGTTAACGAGATCTTTTATGACCTAGTGGCGG  
AAATTAACAGAAAAACCCAGTGCCTGGGAAGGCCCGCAAAAAGTCATCATGTCAGCTGCTT

**ACGCGT**ACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR212319 representing NM\_134346  
Red=Cloning site Green=Tags(s)

MREYKLVVLGSGGVGKSALTVQFVQGI FVEKYDPTIEDSYRKQVEVDAQQCML EILD TAGTEQFTAMRDL  
YMKNGQGFALVYSITAQSTFNDLQDLREQILRVKDTDDVPMILVGNKCDLEDERVVGKEQQNLARQWSN  
CAFLESSAKSKINVNEIFYDLVRQINRKTPVPGKARKKSSCQLL

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online >](#)

## Cloning Scheme:



## Plasmid Map:



ACCN: NM\_134346

ORF Size: 552 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

|                               |   |
|-------------------------------|---|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_134346.3</a></u> , <u><a href="#">NP_599173.2</a></u>   |
| <b>RefSeq Size:</b>           | 555 bp  |
| <b>RefSeq ORF:</b>            | 555 bp  |
| <b>Locus ID:</b>              | 171337  |
| <b>UniProt ID:</b>            | <u><a href="#">Q62636</a></u>   |
| <b>Cytogenetics:</b>          | 7q22  |
| <b>MW:</b>                    | 20.8 kDa  |
| <b>Gene Summary:</b>          | small GTPase; involved in antagonizing the mitogenic and transforming activity of Ras [RGD, Feb 2006]   |