

Product datasheet for **RR212319L3V**

Rap1b (NM_134346) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rap1b (NM_134346) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rap1b
Synonyms:	MGC93206
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_134346
ORF Size:	552 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212319).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134346.3 , NP_599173.2
RefSeq Size:	555 bp
RefSeq ORF:	555 bp
Locus ID:	171337
UniProt ID:	Q62636
Cytogenetics:	7q22
Gene Summary:	small GTPase; involved in antagonizing the mitogenic and transforming activity of Ras [RGD, Feb 2006]



[View online »](#)