

Product datasheet for **RR212340L3V**

Pcdha4 (NM_053933) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pcdha4 (NM_053933) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pcdha4
Synonyms:	Cnr1; cnr5; Pcdha1; Pcdha2; Pcdha5; Pcdha6; Pcdha7; Pcdha8; Pcdha10; Pcdha11; Pcdha12; Pcdha13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053933
ORF Size:	2841 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212340).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053933.1 , NP_446385.1
RefSeq Size:	5241 bp
RefSeq ORF:	2844 bp
Locus ID:	116741
UniProt ID:	Q76718
Cytogenetics:	18p11
Gene Summary:	binds dental epithelial and mesenchymal cells; may play a role in enamel and dentin formation during incisor development [RGD, Feb 2006]



[View online »](#)