

Product datasheet for **RR212348**

Mum111 (NM_001109321) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mum111 (NM_001109321) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mum111
Synonyms:	RGD1565148
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR212348 representing NM_001109321
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGATGGTGAATATGTCCTATGTAGTTGGAAAGACCAATTATGGCCAGCCAAAGTTTTAGATAGGTCTG
 AAAGTCCTTCAGACAGTAAGAGGAAAAAGACATTGTCTGTAGAAGTTGAAATACTTTCACCTGGATGAGAA
 AATTACCGTGGAAAGCAAAGACACAAAAGTCTAACCAAATCTGCTGTCGAAGCCATTATGTCCTCTCTA
 ACAGTGCAGTCAGAGGTCACTTACACCTACGGAGGAACTGCCTATGAAAGATCACTGAAATGGCAC
 TGGAAATCGTGGAAAGAAACAATCAGAGCCACGAAAGCATGGCAGAAGAGCAACATATGGCTACAGC
 ATCTGAAATGTGCCACAGCAGCCGCTGATTCGCCTCCTATAAAAAGTTCCGAAACTCGAAAGCAAC
 CTCCAGGAAGACTCGGCTCCATGTTGCTGTGCTCCGAGAGTGATGATCCCTCTCCGACGAGAAGTTGC
 AGGTGCACACAACCAGTGAGAACATGCCAAGTGAAGTGGAAATGAAGTCAGAAAACCTTAAGCTGCTGCCA
 AACGTACCCTTCTTTTCAGATGATGATGATAAAAAGAAGAAAAGAAAGATTGACATCTCGGCGATC
 ATGCTGTGAATTTATCGCTCAAAGAAGAAAGCGAATTTATTAAGAAGAAAAGTTGTCCCTCATCAG
 AAGATCTCGCTGTACCCAAAGAGGAGTCCCAAGACATCCTCCAGATGCCCGGCTGTTTCTCTGAATG
 CTTTATTGTCTCAGAGAATAACATGGAAGATCCTGGAGAGGGCCCATCAAATCAGAAATCCAGTTTCTAT
 GGCAGCCAAAGTCAGTCTACTGTGGAATCAGATGTAGGTGCTGAGACATCCACTGCAGGGTGCTCAGGGG
 ACTATCAGGTTTCCCTTCTGCCGTGATACAGTCAACAGTGTACTACTCCAGAGACTGGATTTAGA
 AGATCTTGAGGAAGAAGCCGAGCTTCTGGCAAGCTTTTGTCTCTAAATCCTGCCTGTGCGGCTGCGTTA
 GAAAATGATGATGAGGACGACGATGAAGACCTTCCACGTTTATCCTCCGTTATGAGACAGTGCATTTG
 AAACCGCATGATAGTGTGGTTTAAATATCAAAAATACCCATTTTGGCCAGCAGTATCAAAGCATCAG
 GCGGAAAAGAGAGGAAAGCGAGTGTGCTTTTGGTTGAGGCAGACATGAATCCTCAAAGAGAGGCGTTAGA
 GTATCTTTGAGAAGGCTGAAAAATATGACTGTAAAGAGAAACAGGCACTAGTGGAGAAAGCCAGGAGG
 AGTACAAGGAGAGCATTGATTGGTGCCTGCTGATTTGTGACTACCGAGTTAGACTAGGCTGTGGTTC
 TTTTACTGGTTCATTCTTTGAGTATTACGCTGCTGACATCAGTTACCCGGTCAGGAAAATCATCAAACAA
 GATACGTTTCAGAAATATTTTCCAAAGCTATACAATGAAGACGCGGGGAGCAATTGCCTGTGGCATCTC
 ACGCCAAGAGATATCTTCCGAAAATCTCCCTGATCGCATGAAGCCTGCTCGGGACCGAGCTAACAA
 GAACCTTGTAGATTTTATTGTCAATGCAAAGGAACTGAGGACCACCTCCTGGGAATTTTAAAGGGCACA
 AAAAAATCCAAGTGGCTGAAATCATTTCTGAATGCAAAGAGTTTACACCCTGTATTGAAACCTACTTTG
 AAGACGAAGATCAACTGGATGAAGTGGTGAATATCTACAAGAAATTTATAAGCAAATTGACCAGAAGAT
 GCTGACTCTGATAAAAGATGACAAAATTAAGTTTGTCTTGGAAAGTTCTTCTACCCGAAGCAATTTTGC
 TCAATTTCTGCAGTTGATGGCTTAGATTATGAAGCAGCGGAGGCAAAGTATCTAAAGGGACCATCCCTTG
 GCTATAGGGAGAGAGAATTATATGATTCAAAAATCCTATTTGAAAAGAGGCGGAGGTCATTACCAAATGA
 AGGTCGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR212348 representing NM_001109321
Red=Cloning site Green=Tags(s)

```
MDGEYVLCSDKDQLWPAKVLDRSESPDSKRKKTLSVEVEILSLDEKITVESKDTKVLTKSAVEAIMSSL
TVQSEVSFTPEETAYERSLKMALEIVEERTNQSHESMAEEQHMTASENVPQQPPDSPPHKKFRKLESN
LQEDSASMLLCSESDDSLDEKLQVHTTSENMPSEVEMKSENLSGCCQYPSFSDDDDKKEEKKIDISAI
MSVNLSLKEESEFIKEEFVPSSEDLAVPKEESQDILPDAPAVSSECFIVSENNMEDPGEGPSNQNPGSY
GSQSQSTVESDVGAETSTAGCSGDYQVSLPARDTVNSDLLLQRLDLEDLEEEARASGKLLSLNPACAAAL
ENDDEDDDEDLPRFILRYETRAFETGMIYVFKYQKYPFWPAVIKSIIRRKERKASVLLVEADMNPQKRGVR
VSLRRLKKYDCKEKQALVEKAREEYKESIDWCVSLICDYRVRVLCGSGFTGSFFEYAAADISYPVRKIIKQ
DTRFNIFPKLYNEDAGEQLPVASHAKRVSFRKILPDRMKPARDRANKNLVDFIVNAKGTEDHLLGILKGT
KKSKWLKSFNNAKSFTPCIETYFEDEDQLDEVVKYLQEIYKQIDQKMLTIKDDKIKFVLEVLLPEAIIIC
SISAVDGLDYAAEAAYLKGPSLGYRERELYDSKILFEKRRRSLPNEGR
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001109321

ORF Size: 2037 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001109321.1](#), [NP_001102791.1](#)

RefSeq Size: 3796 bp


RefSeq ORF: 2040 bp

Locus ID: 501630

Cytogenetics: Xq32

MW: 77.2 kDa

Product images:



Circular map for RR212348