

Product datasheet for RR212371

Kctd18 (NM_001106914) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Kctd18 (NM_001106914) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Kctd18
Synonyms: RGD1564820
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR212371 representing NM_001106914
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCTGGCCACGAAGCAGAGGAGATGCGAGATATCCTTCGACTGAACGTGGGCGGCTGCATTTACACAG
 CTCGGCGAGAGTCTTTGTGCCGCTTAAGGACTCTATGCTGGCATCTATGTTAGCGGTCTGTTCCCTCT
 CAAAACAGATGAATCTGGGGCGTGTATTATTGACCGTGATGGACATCTATTCAAATACATTTGGACTAC
 CTCCATGGGGAGGTTAGACCCCTCGGACGAGCAGACCCGAACCGCCCTGCAGGAAGAGGCTGACTACT
 TTGGCATCCCTTACCCATACAGCCTGTCTGACCACTTGGCCAATGAGATGGAGACATATTCTTTAAGGTC
 AAATATAGAACTTAAAAAGGCTTTAACAGACTTCTGTGATTCTATGGCTTAGTTTGAATAAACCAACA
 GTTTGGGTCTCCACTATCTTAACACATCTGGTGAAGCTGTGAAAGTAGAATTATTGGCGTCTATTCTA
 CAAAACTGATGGGACAGATGCCATTGATAAGCAGCTGGGAGGAAGAATTCATAGTAAAAGCATTTCAA
 AAGGGAGGCGGGGAACAATGTCCAGTATATTGGAGCTACTACTCAGTAGCCGAGTTGAAGAAAAATGATG
 GACGCCCTTTGATGCTTGGGAAGGAAAAGGTGAAGCTACTGGCGTGTGCCTCATGAGCTGATAGAGTGTT
 GGACTTTGGAAGAACGGCCTTTACATGGAAGCCTGCGTCATATGGCCCCAATTCGAAAGAGGCGGCTGAT
 AGCCTTTGGGGAAGAGGATGAAGGTGTGAAGTGTAAAAGTGGTCCCTAAGCCAGTCAGATTTTGGGCCCT
 TCTACAAGTACCCAAATTAAGTCAAGAACTCTGCTTCAGTGAGGGTGTACCTGCCAATGCTGCCCAGA
 TCTCCTCGCAGGCTACAGCAAACCGCTCTCAAAGAGGCAGTTGTGAAAGGCAGCCAGTCTGTGGC
 CACAGGGGCATCTGGACATGCTCCTGCTTCTCCTAATGCCAGAGCGCTGAAAACCAAGCCACACAGCCA
 CCTCCTGCCAAGGTGCTCCTCTGTGATAAGAAGTCTCCACCCATCGAGTGATCAAGCTGAGGAGGACTC
 CACTGTGTGTGGGCCTTCCTGCCAACCCCTCAGCAGGGAGCCAGGCCAGCTTCCCCAGCAGCCCTC
 TCCAGAAGCTTCCAACACTTTGGGCGTGCAGACTGAGAACAGGAAAGACCAGCCTGATGGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR212371 representing NM_001106914
Red=Cloning site Green=Tags(s)


MAGHEAEEMRDILRLNVGGCIYTARRESLCRFKDSMLASMFSGRFPLKTDESGACIIDRDGHLFKYILDY
 LHGEVQTPSDEQTRTALQEEADYFGIPYPYSLSDHLANEMETYSLSNIELKKALTDFCDSYGLVCNKPT
 VWVLHLYLNTSGASCESRIIGVYSTKTDGTAIDKQLGGRIHKSIFKREAGNNVQYIWSYYSVAELKMM
 DAFDAWEGKGVSYWRVPHELIECWTLERPLHGSLRHMAPIRKRRLIAFGEEDEGVNCKTGPKPVRFLGP
 STSTQIKVKNSASVRVSPANAAQISSQATANRSQRGSCGKAAQCSVATGASGHAPASNAQSAENQATQP
 PPAKVLLSDKKSPPHRVIKLRRRTPLCVGPSPPTPSAGSQASFPQQPSPEASNTLVGQVTENRKDQPDG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106914

ORF Size:	1251 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001106914.1 , NP_001100384.1
RefSeq Size:	5901 bp
RefSeq ORF:	1254 bp
Locus ID:	301436
Cytogenetics:	9q31
MW:	45.9 kDa