

Product datasheet for **RR212389L3V**

Mif (NM_031051) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mif (NM_031051) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Mif
Synonyms:	MGC72801
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031051
ORF Size:	345 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212389).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031051.1 , NP_112313.1
RefSeq Size:	551 bp
RefSeq ORF:	348 bp
Locus ID:	81683
UniProt ID:	P30904
Cytogenetics:	20p12
Gene Summary:	inhibits random migration of macrophages and is involved in the pathogenesis of several inflammatory diseases [RGD, Feb 2006]



[View online »](#)