

## Product datasheet for **RR212412**

### **Mbl2 (NM\_022704) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Mbl2 (NM\_022704) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Mbl2  
**Synonyms:** Ab2-001; Ab2-011  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR212412 representing NM\_022704  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCCCTGTTACATCCTTCCTTCTGCTCTGCGTGCTCACGGCAGTCTATGCCGAGACCTTAACCGAAG  
GGGCTCAAAGTAGCTGCCCTGTGATTGCCTGCAGTTCTCCGGCCTGAACGGCTTCCAGGCAAAGATGG  
ACACGACGGTGCCAAGGGAGAAAAGGGAGAACCGGGTCAAGGCCTCAGAGGCTTGCAGGGCCCTCCTGGA  
AAAGTAGGACCTGCAGGGCCCCAGGGAATCCTGGGTCAAAGGAGCAACGGGACCAAAGGAGACCGTG  
GAGAGAGTGTAGAATTTGATACTACCAACATTGATTTAGAAATTGCAGCCCTGCGATCGGAGCTGAGAGC  
TATGAGAAAGTGGGTGCTCCTTTCTATGAGTGAATGTTGAAAGAAGTACTTCATGAGCAGTGTAGAG  
AGGATGCCCTTAACAGAGCGAAGGCTCTGTGCTCCGAACTCCAGGGCACTGTGGCCACTCCCAGGAATG  
CTGAGGAAAATAGGGCCATCCAGAATGTGGCAAAGATGTTGCCTTCTTGGGCATAACGGACCAGAGGAC  
TGAAAACGTTTTTGGAGACCTGACAGGAAACAGAGTGCCTACACTAACTGGAATGAGGGTGAGCCCAAC  
AATGTGGGCTCTGGGAAAACGTGTGGTGTCTTGACAAATGGGAAGTGAATGACGTTCTTGTCTGTG  
ATTCTTTTTGGTAGTTTGTGAATCTCTGAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR212412 representing NM\_022704  
Red=Cloning site Green=Tags(s)

MSLFTSFLLLCVLTAVYAETLTEGAQSSCPVIACSSPGLNGFPGKDGHDGAKGEKGEPEGQLRGLQGPPG  
 KVGPPAGPPGNPGSKGATGPKGDRGESVEFDTTNIDLEIAALRSELRAMRKWVLLSMSENVGKKYFMSSVR  
 RMPLNRAKALCSELQGTVATPRNAEENRAIQNVAKDVAFLGITDQRTENVFEDLTGNRVRVYTNWNEGEPN  
 NVGSGENCVVLLTNGKWNDVPCSDSFLVVFCEFS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_022704  
**ORF Size:** 732 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_022704.2</a> , <a href="#">NP_073195.2</a>
<b>RefSeq Size:</b>	874 bp
<b>RefSeq ORF:</b>	735 bp
<b>Locus ID:</b>	64668
<b>UniProt ID:</b>	<a href="#">P08661</a>
<b>Cytogenetics:</b>	1q52
<b>MW:</b>	26 kDa
<b>Gene Summary:</b>	a lectin that binds mannose and N-acetylglucosamine [RGD, Feb 2006]