

## Product datasheet for **RR212413L4V**

### Hmgn3 (NM\_001007020) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hmgn3 (NM_001007020) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hmgn3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001007020
ORF Size:	285 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212413).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001007020.2</a> , <a href="#">NP_001007021.1</a>
RefSeq Size:	1434 bp
RefSeq ORF:	288 bp
Locus ID:	113990
UniProt ID:	<a href="#">Q66H40</a>
Cytogenetics:	8q31



[View online »](#)

**Gene Summary:**

Binds to nucleosomes, regulating chromatin structure and consequently, chromatin-dependent processes such as transcription, DNA replication and DNA repair. Affects both insulin and glucagon levels and modulates the expression of pancreatic genes involved in insulin secretion. Regulates the expression of the glucose transporter SLC2A2 by binding specifically to its promoter region and recruiting PDX1 and additional transcription factors. Regulates the expression of SLC6A9, a glycine transporter which regulates the glycine concentration in synaptic junctions in the central nervous system, by binding to its transcription start site. May play a role in ocular development and astrocyte function (By similarity).[UniProtKB/Swiss-Prot Function]