

Product datasheet for RR212437

Pde7b (NM_080894) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pde7b (NM_080894) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pde7b
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR212437 representing NM_080894 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTCTTGTTAATGGTTGAGAGGTGTGGCGAAGTCTTGTGTTGAGAATCCTGAGCAGAATGTCAAATGTG
TTTGCATGCTAGGAGATGTACGACTAAGGGTCTAGCAGGGGTTCTGCCGAACGCCGGCTCCTACCC
GTTTCATTGACTTCGGTCTACTTAAACAATAACAACACTCAGGGGAAATTGGCAGCAAGAAAAAGGTGAGG
AGACTGTTAAGTTTCCAAGGCACTTCCATGAATCTAGGCTGCTCCGGGGATGACACCGCAGGCCCC
TCCACCTGCTGGACGAAGACTACCTTGGACAAGCAAGGCATATGCTCTCCAAAGTTGGAATGTGGGACTT
TGACATTTTCTGTTGATCGCTTGACAAATGGGAACAGTCTGGTAACTCTGTTGTGTCACCTCTTCAAC
TCCCATGGACTCATCCACCATTCAAGCTTGACATGGTGACCTTACACAGGTTTTGGTTATGGTTCAGG
AAGATTACCACGGCCACAACCCGTACCACAATGCTGTTTCATGCAGCTGACGTACCCAGGCCATGCACTG
TACTTGAAGGAGCCAAAGTTGGCGAGCTTCTCACACCTCTGGACATCATGCTTGGACTATTGGCTGCA
GCAGCTCATGACGTGGACCACCCAGGGTGAACCAGCCATTTTGTGATCAAACTAACCACCCTTGCCA
ACTTGTATCAGAACATGTCGGTACTGGAGAATCACCCTGGCGGTCTACAATCGGCATGCTTCGGGAATC
AAGGCTCCTTGCTCATTGCCAAAGGAAATGACACAGGATATCGAACAGCAGCTGGGCTCCTGATCTTG
GCCACGGACATCAACAGACAGAATGAGTTCCTGACCCGCTTAAAAGCTCACCTCCACAATAAGGATTTGA
GACTGGAAAATATACAGGACAGACACTTATGCTTCAAGTGCAGTGCCTTGAAGTGCAGTGCATTTGCAATCC
TTGTGCAATCTGGGAGATGAGCAAGCAGTGGAGTGAAGAGTCTGCGAAGAATTCTACAGGCAAGGTGAC
CTTGAACAGAAGTTTGAAGTGAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGT
AAATTGGTTTCATGACTTACATCGTGGAGCCACTGTTCCGGGAGTGGCCCGGTTCACTGGGAACAGCAC
CCTGTCCGAGAGCATGCTAAACCATCTCGGCACAACAAGCCAGTGAAGAGCCTGCTGTCCAATCAG
CACAGACGAGGGCAGCGGCCAGGACCCAGCGGCACGGCGCTGAGACCCTGGAGCAGACAGAAGGCC
CCACGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)


Protein Sequence: >RR212437 representing NM_080894
Red=Cloning site Green=Tags(s)

MSCLMVERCGEVLFFENPEQNVKVCMLGDVRLRGQTGVAERRGSYPFIDFRLLNNTTHSGEIGSKKKVR
 RLLSFQRHFHESRLLRGMTPQAPLHLLDEDYLQARHMLSKVGMWDFDIFLFDRLTNGNSLVTLCHLFN
 SHGLIHHFKLDMVTLHRFLVMVQEDYHGHNPYHNAVHAADVTQAMHCYLKEPKLASFLTPLDIMLGLLAA
 AAHDVDHPGVNQPFLLIKTNHHLANLYQNMVLENHHWRSTIGMLRESRLLAHLPKEMTQDIEQQLGSLIL
 ATDINRQNEFLTRLKAHLHNKDLRLENIQDRHFMLQIALKCADICNPCR IWEMSKQWSEVCEEFYRQGD
 LEQKFELEISPLCNQKDSIPSIQIGFTYIIVEPLFREWARFTGNSTLSESMNLHLAHNAQWKSLLSNQ
 HRRRSGSQDPAGTAPETLEQTEGATP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_080894
ORF Size:	1338 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_080894.2 , NP_543170.1
RefSeq Size:	1753 bp
RefSeq ORF:	1341 bp
Locus ID:	140929
Cytogenetics:	1p12
MW:	51.5 kDa
Gene Summary:	catalyzes the conversion of adenosine 3',5'-cyclic phosphate to adenosine 5'-phosphate [RGD, Feb 2006]