

## Product datasheet for **RR212472L3V**

### Slc26a3 (NM\_053755) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc26a3 (NM_053755) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc26a3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053755
ORF Size:	2271 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212472).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053755.2</a> , <a href="#">NP_446207.1</a>
RefSeq Size:	2850 bp
RefSeq ORF:	2274 bp
Locus ID:	114629
UniProt ID:	<a href="#">Q924C9</a>
Cytogenetics:	6q16
Gene Summary:	transmembrane chloride/bicarbonate exchanger; plays a role in mediating electrolyte and fluid absorption in the colon [RGD, Feb 2006]



[View online »](#)