

## Product datasheet for RR212478L4

### Dgkz (NM\_031143) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Dgkz (NM_031143) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Dgkz
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212478).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



RefSeq ORF: 2790 bp  
Locus ID: 81821  
UniProt ID: [O08560](#)  
Cytogenetics: 3q24  
Gene Summary: catalyzes the conversion of diacylglycerol to phosphatidic acid [RGD, Feb 2006]