

Product datasheet for RR212482L4V

Chmp3 (NM_172331) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Chmp3 (NM_172331) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Chmp3
Synonyms:	Vps24
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
ACCN:	NM_172331
ORF Size:	669 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212482).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172331.1 , NP_758834.1
RefSeq Size:	831 bp
RefSeq ORF:	672 bp
Locus ID:	282834
Cytogenetics:	4q31
Gene Summary:	binds phosphatidylinositol 3,5- and phosphatidylinositol-4,5-bisphosphates; involved in endosome homeostasis [RGD, Feb 2006]



[View online »](#)