

Product datasheet for RR212492

Efna4 (NM_001107692) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Efna4 (NM_001107692) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Efna4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR212492 representing NM_001107692
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCGGCTGTTGCCCTGCTGCGGACTGTGCTCTGGGCCGCGCTGCTCGGCTCGCGCCTGCGGGGTGCT
 CCAGCCTCCGCACTCTACTGGAAGTCCACTAAGCCAGGTTGCTTCGAGGAGATGCCGTGGTGA
 GCTGGGCTCAACGATTACCTAGACATCTTCTGCCACATTATGAGAGCCAGGGCCCCAGAGGGCCCC
 GAAACGTTTGCATTATACATAGTGGACTGGTCAAGGCTGCAAGGCAGAGGGGCAAAATGCCT
 TCCAGCGCTGGAATTGCACGCTACCTTTTGGCCCTTTTGTCCCTGTTTCGATTCTCCGAAAAGATTACAGCG
 CTTACACCGTTCCCGCTGGGCTTCGAGTCTTGCCTGGAGAGACTTACTACTACATCTCGGTGCCAACT
 CCGGAGAGTCTGGCCAGTGCCTGAGACTCCAGGTGTCTGTCTGCTGCAAGGAGGACAGATCCGAGTCAG
 CTCATCTGTCCGGAGCCCTGGAGAAAGTGGCACATCTGGGTGGCGGGGAGGACACGCCCCAGCCCTCT
 GTCTCTCTGCTGTTGCTGCTGCTCCCAATACTCCGTCTTCTGCGAGTTCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR212492 representing NM_001107692
 Red=Cloning site Green=Tags(s)

MRLLPLLRTVLWAALLGSRLRGCSLRSIYWNSTNPRLLRGDAVVELGLNDYLDIFCPHYESPGPPEGP
 ETFALYIVDWSGYEACKAEGANAFQRWNCTLPFAPFVVPVRFSEKIQRFTPFPLGFELPGETYYYYISVPT
 PESPGQLRLQVSVCKEDRSESAHPVGSPEGSTSGWRGGHAPSPLSLLLLLLLLPILRLLRVL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001107692

ORF Size: 612 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001107692.1, NP_001101162.1
 RefSeq Size: 1199 bp
 RefSeq ORF: 615 bp
 Locus ID: 310643
 Cytogenetics: 2q34
 MW: 22.8 kDa

Product images:



Circular map for RR212492