

Product datasheet for **RR212516**

Dhx36 (NM_001107678) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dhx36 (NM_001107678) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dhx36
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR212516 representing NM_001107678
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAGCTACGACTATCATCAGAGCTGGAGCCGCGACGGGGGCCACGGGGCTCCGGCCAGGGCTCTGGCG
 GCGGCGGTGGAGGAGCCGGGGCTCCGGCGGCGCGAGGGGGCCGCGGGCGGGCCGGCATCCGGC
 CCACCTTAAGGGTCGCGAGATAGGCCTATGGTACGCAAAAAAGCAAACGAGAAGAACAAGGAGGCTGAG
 AGACAAGAGAGAGCTGTAGTGCACATGGATGAACGGCGAGAAGAGCAAATTGTGCAGCTGCTAAATTCAG
 TCCAAGCAAAGAATGATAAAGATTCAGAAGCACAGATATCCTGGTTTGTCTCTGAGGATCATGGGTATGG
 TACTGAAGTTTCTTCAGAGAAAAAATAAACTCAGAGAAGAACTTGACAACCAGGAAAAGAAATTTGCTG
 AACCAAGAAAAAAGACATATAGGATCACAGACAAATCATATATTGACCGAGATTCTGAGATTTTATTGC
 AACAAAATGAACCAACCTAGGCTTAGACCAACAATTATTGGAAGATTTACAAAAGAAAAAACTGACCC
 TCGGTATATAGAGATGCAGCGTTTCAGAAAAAGCTGCCTTCATATGGAATGCAGAAGGAGCTGGTAAAT
 TTAATCAATAACCATCAGGTGACAGTGATAAGTGGTGAGACTGGTTGTGGCAAAACCACTCAGGTTACTC
 AGTTCATCTTGGATAACTACATCGAAAGAGGAATAGGGTCTGCCTGCAGAATAGTGTGACTCAGCCAAG
 AAGAATTAGTGCCATTTAGTTGCTGAGAGGGTGGCTGCAGAAAGGGCAGAGTCTTGTGGCAATGGTAAT
 AGTACTGGATACCAGATTCGTCTCCAAAGTCGGTTGCCAAGGAAACAAGGCTCTATCCTATACTGCACAA
 CAGGAATCATTCTTCAGTGGCTCCAGTCAGACTCACGTTTGTCCAGTGTTAGTCAATTTGACTTGTGA
 AATTCATGAAAGGAATCTACAGTCAGATGTTTTAATGACTGTTATTAAGATCTCCTCCATTTTCGATCG
 GATCTCAAAGTAATATTGATGAGTGAACCTTGAATGCTGAGAAATTTTCAGAATATTTTGGTAACTGTC
 CAATGATACATATACCTGGGTTACTTTCCAGTTGTGGAATATCTTTTGGAAAGATATCTGAAAAAAT
 AAGATATTTTCCAGAACAAAAAGAACATAGATCCAGTTCAAAGGGGTTTCATGCAGGGTCATGTAAT
 AGACAAGAAAAAGAAGAAAAAGAGGCGATCTATAAAGAACGCTGGCCAGCTTATATAAGGAACTACAGA
 CAAGATACTCTGCAAGTACCATAGATGTTTTGAAATGATGGATGATGATAAAGTTGATCTGAATTTGAT
 TGCTGCCCTTATTCGATACATTGTTTTGGAAGAAGAGGACGGTGCATATTGGTCTTTCTACCAGGCTGG
 GACAATATCAGTACTTTGCATGATCTCTTGATGTCACAAGTGATGTTAAATCAGATAGGTTTCTTATTA
 TACCTTTACATTCAGTACTGATGCCACTGTAACCAGACACAGGATTTAAAAAACTCCTCCTGGTGTTCG
 GAAAATAGTAATTGCTACCAACATTGCAGAGACTAGCATACCATAGATGATGTGGTCTATGTAATAGAT
 GGAGGGAAAATTAAGAAACACACTTTGACACACAGAACAACATCAGTACCATGTCTGCTGAGTGGGTTA
 GTAAGCTAATGCCAAACAGAGGAAAGGAGAGCTGGAAGAGTTACGCTGGTCAATGTTATCATCTGTA
 TAATGGTCTTAGAGCAAGTCTTCTAGATGACTACCAACTACCAGAAATTTTGAGAACTCCTTTGGAAGAA
 CTTTGTGTTGCAAATAAAGATCTTGAGGCTGGGTGGAATTGCTTATTTTCTGAGTAGGTTAATGGATCCAC
 CATCAGATGAGGCAGTAGTGCTCTCCAATAAACATCTAATGGAAGTGAAGTCTTTGGATAAGCAAGAAGA
 ATTGACACCTCTTGAGTCCATTTAGCCCGACTGCCTGTTGAGCCACATATTGAAAAATGATTCTCTTT
 GGAGCTTTGTTCTGTTGCTTAGATCCAGTCCCTACCATTGCTGCCAGTCTCAGCTTTAAAGATCCCTTTG
 TCATTCATTGGGAAAAGAAAAGATTGCAGATGCAAGAAGAAAGGAACTGGCAAAGGAACTAGAAGTGA
 TCATCTGACAGTTGTAATGCGTTTGGGGCTGGGAAGAAGCTAAACGGCGTGGTTTCAGATATGAAAAA
 GACTACTGCTGGGAATATTTCTCTCTTCAAACACACTACAGATGCTGCATAACATGAAGGGACAGTTTG
 CTGAGCATCTTCTGGAGCTGGATTTGTAAAGCAGTAGAAGTCCCAAAGATCCAAAAGCTAATATAAATTC
 AGATAATGAGAAGATCATTAAAGCTGTCATCTGTGCTGGTTTATATCCCAAAGTTGCTAAAAATCCGACTA
 AATTTGGGTAAAAAAGGAAAAATGGTGAAGTTTACACAAAGTCTGATGGTCTGGTTTCTATTATCCTA
 AGTCTGTCAATGTGGAGCAGACAGACTTCCACTACAAGTGGCTTATCTATCACCTGAAGATGAGAACAAG
 CAGTATACTTGTACGACTGCACAGAAGTCTCACCATACTGCCTCTTGTCTTTGGAGGAGATATTTCC
 ATCCAGAAAAGATAAGGATCAGGAAATATTGCTGTAGACGAATGGATTGTGTTTCAGTCTCCAGAAAAGAA
 TTGCTCATCTTGTAAAGGACTAAGAAAGGAACTGGATATCCTTCTACAAGAGAAGATTGAGTGCCTCA
 TCCTGTAGATTGGAACGACACCAATCCAGAGACTGTGCAAGTACTGTACGGATTCTAGACTTGATCAAA
 ACACAAGAAAAGGCTATTTCAAGGAACTTGCCACCACGATCACAGGATGGATATTACAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR212516 representing NM_001107678
 Red=Cloning site Green=Tags(s)


MSYDYHQSWSRDGGPRGSGQSGGGGGSRGSGGGGGGRGRGRHPAHLKGREIGLWYAKQTQKNKEAE
 RQERAVVHMDERREEQIVQLLNSVQAKNDKDSEAIQSWFAPEDHGYGTEVSSEKKINSEKLDNQEKKLL
 NQEKKTYRITDKSYIDRSEYLLQQNEPNLGLDQQLLEDLQKKTDPRYIEMQRFRKKLPSYGMQKELVN
 LINNHQVTVISGETGCGKTTQVTQFILDNYIERGIGSACRIVCTQPRRISAISSVAERVAEAERESCGNGN
 STGYQIRLQSRLPRKQGSILYCTTGIIQLWLQSDSRLSSVSHIVLDEIHERNLQSDVLMTVIKDLLHFRS
 DLKVILMSATLNAEFSEYFGNCPMIHIPGFTFPVVEYLLEDIIEKIRYFPEQKEHRSQFKRGMQGHVN
 RQEKEEKEAIYKERWPAYIKELQTRYASTIDVLEMMDDDKVDLNLIAALIRYIVLEEDGAILVFLPGW
 DNISTLHDLLMSQVMFKSDRFLIIPHLMLPTVNQTQVFKKTPPGVRKIVIATNIAETSITIDVVVVID
 GGIKIKETHFDTQNNISTMSAEWVSKANAKQRKGRAGRVPGHCHLYNGLRASLLDDYQLPEILRTPLEE
 LCLQIKILRLGGIAYFLSRLMDPPSDEAVVLSIKHLMEL SALDKQEEL TPLGVHLARLPVEPHIGKMILF
 GALFCCLDPVLTIAASL SFKDPFVIPLGKEKIADARRKELAKETRSDHL TVVNAFEGWEEAKRRGFRIEK
 DYCWEIFLSSNTLQMLHNMKGQFAEHL L GAGFVSSRSPKDPKANINSDNEKIIKAVICAGLYPKVAKIRL
 NLGKKRKMVKVHTKSDGLVSIHPKSVNVEQTDHFYNWLIYHLKMRSSIYLYDCTEVSPYCLLFFGGDIS
 IQKDKDQEIIAVDEWIVFQSPERIAHLVKGLRKELDILLQEKIECPHPVDWNTKSRDCAVLSAILDLIK
 TQEKAI PRNLPPRSQDGYYS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001107678

ORF Size: 3000 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001107678.1](#), [NP_001101148.1](#)

RefSeq Size: 5664 bp

RefSeq ORF: 3003 bp

Locus ID: 310461

UniProt ID: [D4A2Z8](#)

Cytogenetics: 2q31

MW: 113.8 kDa

Gene Summary: Multifunctional ATP-dependent helicase that unwinds G-quadruplex (G4) structures (By similarity). Plays a role in many biological processes such as genomic integrity, gene expression regulations and as a sensor to initiate antiviral responses (PubMed:23651854). G4 structures correspond to helical structures containing guanine tetrads (By similarity). Binds with high affinity to and unwinds G4 structures that are formed in nucleic acids (G4-ADN and G4-RNA) (By similarity). Plays a role in genomic integrity. Converts the G4-RNA structure present in telomerase RNA template component (TREC) into a double-stranded RNA to promote P1 helix formation that acts as a template boundary ensuring accurate reverse transcription (By similarity). Plays a role in transcriptional regulation. Resolves G4-DNA structures in promoters of genes, such as YY1, KIT/c-kit and ALPL and positively regulates their expression (By similarity). Plays a role in post-transcriptional regulation. Unwinds a G4-RNA structure located in the 3' UTR polyadenylation site of the pre-mRNA TP53 and stimulates TP53 pre-mRNA 3'-end processing in response to ultraviolet (UV)-induced DNA damage (By similarity). Binds to the precursor-microRNA-134 (pre-miR-134) terminal loop and regulates its transport into the synapto-dendritic compartment (PubMed:23651854). Involved in the pre-miR-134-dependent inhibition of target gene expression and the control of dendritic spine size (PubMed:23651854). Plays a role in the regulation of cytoplasmic mRNA translation and mRNA stability. Binds to both G4-RNA structures and alternative non-quadruplex-forming sequence within the 3' UTR of the PITX1 mRNA regulating negatively PITX1 protein expression. Binds to both G4-RNA structure in the 5'-UTR and AU-rich elements (AREs) localized in the 3' UTR of NKX2-5 mRNA to either stimulate protein translation or induce mRNA decay in an ELAVL1-dependent manner, respectively. Binds also to ARE sequences present in several mRNAs mediating exosome-mediated 3'-5' mRNA degradation. Involved in cytoplasmic urokinase-type plasminogen activator (uPA) mRNA decay (By similarity). Component of a multi-helicase-TICAM1 complex that acts as a cytoplasmic sensor of viral double-stranded RNA (dsRNA) and plays a role in the activation of a cascade of antiviral responses including the induction of proinflammatory cytokines via the adapter molecule TICAM1. Required for the early embryonic development and hematopoiesis. Involved in the regulation of cardioblast differentiation and proliferation during heart development. Involved in spermatogonia differentiation. May play a role in ossification (By similarity).[UniProtKB/Swiss-Prot Function]