

Product datasheet for **RR212522L3V**

Bglap (NM_013414) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bglap (NM_013414) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Bglap
Synonyms:	Bglap2; Bgp; Bgpr; Bgpra
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013414
ORF Size:	297 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212522).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013414.1 , NP_038200.1
RefSeq Size:	487 bp
RefSeq ORF:	300 bp
Locus ID:	25295
UniProt ID:	P04640
Cytogenetics:	2q34
Gene Summary:	highly conserved protein associated with mineralized bone matrix; protein is secreted by calcified tissues and is regulated by vitamin D3 [RGD, Feb 2006]



[View online »](#)