

## Product datasheet for **RR212553L3V**

### Park2 (Prkn) (NM\_020093) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Park2 (Prkn) (NM_020093) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Prkn
Synonyms:	Park; Park2; Prkn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020093
ORF Size:	1395 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212553).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020093.1</a>
RefSeq Size:	1564 bp
RefSeq ORF:	1398 bp
Locus ID:	56816
UniProt ID:	<a href="#">Q9JK66</a>
Cytogenetics:	1q11
Gene Summary:	localized to the cytoplasm in neurons of many regions of the brain; has N-terminal ubiquitin and ring-IBR-ring finger domains [RGD, Feb 2006]



[View online »](#)