

Product datasheet for **RR212594**

Ddx3y (NM_001108858) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ddx3y (NM_001108858) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ddx3y
Synonyms:	RGD1309586
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RR212594 representing NM_001108858
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGATCGCC

ATGAGTCATGTGGCAGAGGAAGATGAGCTTGGGCTGGACCAGCAGTTGGCTGGCCTAGACCTGACGTCTC
CAGACAGCCAGAGCGGGCAAGCACAGCCAGCAAAGGGCGCTACATCCCTCCTCACCTGCGGAACAGGGA
GGCTACGAAAGCGTCCTACGACAAAGACAGCTCACGGTGGAGTAAAGACAAGGATGCCTACAGCAGCTTT
GGTTCCCGCAGTGACACGAGAGCAAGTCGAGTTTCTTCAGCGACCGGGAAGTGCTCCAGAGGAAGGT
TTGATGACCGGGCGGAGTGACTATGAGGGCGTTGGCAGCCGTGGAGGTAGAAGTGGCTTTGGCAAATT
TGAGCGGGGCGGGAACAGTCGCTGGTGTGACAAAGCAGACGAGGATGACTGGTCTAAACCACTCCCACCA
AGCGAACGCTTGAACAGGAGCTCTTTTCTGGAGGCAATACTGGGATTAACCTTGAGAAAATATGATGACA
TTCCAGTTGAAGCGACAGGGAATAACTGCCCTCCACATATTGAAAGTTTCAGCGACGTCGAGATGGGAGA
AATTATTATGGGAAACATCGAGCTCACTCGTTATACACGCCAACTCCGGTACAAAAGCATGCTATTCTCT
ATTATCAAGGAGAAAAGAGACTTGATGGCTTGTGCTCAAACAGGGTCTGGAAAAACGGCAGCTTTTCTTC
TACCTATCTTGAGTCAGATCTATACAGATGGTCCAGGTGAGGCTTTGAGGGCCATGAAGGAGAATGGAAA
GTATGGGCGTCGCAACAGTATCCGATCTCCTTGGTTTTGGCTCCGACCAGAGAATTGGCAGTGCAGATC
TATGAGGAGGCCAGGAAATCTCATACCGATCTCGAGTTCGACCTTGTGTGGTTTATGGTGGAGCTGATA
TTGGGCAGCAGATTTCGAGACTTGAACGTGGCTGTCACTTGTAGTAGCTACTCCAGGACGGCTAGTGGA
CATGATGGAGAGGAAAGATCGGACTGGACTTCTGCAAATACTGGTGTAGACGAAGCGGATCGAATG
TTGGACATGGGCTTTGAGCCTCAGATACGAAGAATAGTTGAACAAGACACCATGCCTCCAAAAGGTGCC
GTCACACTATGATGTTTAGCGCTACTTTTCCTAAGGAAATACAGATGCTGGCTCGGATTTCTTGGATGA
ATATATCTTTCTGGCCGTAGGAAGAGTCGGATCTACTTCTGAGAACATCACACAGAAAAGTTGTGGGGTG
GAAGAAGCAGACAAGCGGTCATTTCTGCTTACCTCCTGAATGCGACAGGCAAGGATTCCTGATATTGG
TGTTCTGGAAACAAAAAGGGTGCAGATTCTCTGGAGATTCTTATACCATGAAGGATATGCTTGTAC
CAGTATTCATGGAGACCGATCTCAGCGAGATAGGGAAGAGGCCCTTACCAATTCGATCAGGAAAAAGT
CCAATTCTCGTGGCTACAGCAGTAGCAGCAAGAGGACTGGATATTTCCAATGTGAAACATGTCATCAATT
TCGACTTGCTAGTGACATCGAAGAATATGTGCATCGAATTGGCCGTACAGGTCGTGTTGGAAACCTTGG
CCTTGCCACCTCATTCTCAACGAGAGGAACATAAATATTACCAAGATTTATTGGATCTCTTGTGAA
GCAAAACAAGAAGTGCCATCTTGGCTAGAAAACATGGCTTTTGAACATCATTACAAGGGTGGCAGTCGTG
GACGGTCCAAGAGCAGATTTAGTGGAGGATTTGGTGTAGAGACTATCGACAAAGTAGCGGTGCCAGCAG
TTCCAGCTTACAGAGTGGCCGTGCTAGCAACAGTCGACAGTGGTGGAGGCAGCCATGGCAGCAGCAGAGGG
TTTGGTGGAGGTAGCTATGGAGGCTTTACAATAGTGTGATGGATATGGAGGAAATTATAGCTCCCAGGGGG
TTGACTGGTGGGTAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR212594 representing NM_001108858
Red=Cloning site Green=Tags(s)

```
MSHVAEEDELGLDQQLAGLDLTSPDSQSGASTASKGRYIPPHLRNREATKASYDKDSSRWSKDKDAYSSF
GSRSDTRAKSSFFSDRGSGSRGRFDDRGRSDYEGVGSRRGSRGFKFERGGNSRWCDKADEDDWSKPLPP
SERLEQELFSGGNTGINFEKYDDIPVEATGNCCPPHIESFSDVEMGEIIMGNIELTRYTRPTPVQKHAIP
IIKEKRDLMACAQTGSGKTAFFLLPILSQIYTDGPGREALRAMKENGKYGRKQYPIISLVLAPTRELAVQI
YEEARKFSYRSRVRPCVVYGGADIGQQIRDLERGCHLLVATPGRLVDMMERGKIGLDFCKYLVLDEADRM
LDMGFEPQIRRIVEQDTMPPKGVRRHTMMFSATFPKEIQMLARDFLDEYIFLAVGRVGTSENITQKVVVV
EEADKRSFLLDLLNATGKDSLILVFVETKKGADSLDFLYHEGYACTSIHGDRSQRDREEALHQFRSGKS
PILVATAVAARGLDISNVKHVINFDLPSEIEEYVHRIGRTGRVGNLGLATSFNERNINITKDLLLLEVE
AKQEVPSWLENMAFEHHYKGGSRGRSKSRFSGGFARDYRQSSGASSSSFSSGRASNSRSGGSHGSSRG
FGGGSYGGFYNSDGYGGNYSSQGVDDWGN
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001108858

ORF Size: 1977 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001108858.1](#), [NP_001102328.1](#)

RefSeq Size: 3170 bp


RefSeq ORF: 1980 bp

Locus ID: 364073

Cytogenetics: 13q26

MW: 73 kDa

Product images:



Circular map for RR212594