

Product datasheet for **RR212609**

Bmp4 (NM_012827) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bmp4 (NM_012827) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Bmp4
Synonyms: BOMPR4A
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR212609 representing NM_012827
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGATTCCTGGTAACCGAATGCTGATGGTCGTTTTATTATGCCAAGTCCTGCTAGGAGGCGCGAGCCATG
CTAGTTTGATACCTGAGACCGGGAAGAAAAAGTCGCCGAGATTCAGGGCCACCGGGAGGACGCCGCTC
AGGGCAGAGCCATGAGCTCCTGCGGGACTTCGAGGCGACACTTCTGCAGATGTTTGGGCTCGCCCGCT
CCGCAGCCGAGCAAGAGCGCCGTATCCCGGATTACATGAGGGATCTTTACCGGCTCCAGTCTGGGGAGG
AGGAGGAAGAAGAGCAGAGCCAGGGAACCGGGCTTGTAGTACCCTGAGCGTCTGCCAGCAGAGCCAAAC
TGTGAGGAGTTTCCATCACGAAGAATCTGGAGAATCCCAGGGACCAGTGAGAGCTCTGCTTTTCGT
TTCTTCTTTAACCTCAGCAGCATCCAGAGAATGAGGTGATCTCCTCTGCAGAGCTCCGGCTATTTCCGG
AGCAGGTGGACCAGGGCCCTGACTGGGAACAGGGCTTCCACCGTATAAACATTTATGAGGTTATGAAGCC
CCCAGCAGAAATGGTGCCTGGACACCTCATCACGACTACTGGACACCAGACTAGTCCATCACAATGTG
ACACGGTGGGAACTTTTCGATGTGAGCCCTGCAGTCTTCGATGGACCCGGGAAAAGCAACCAACTATG
GGCTGGCCATTGAGGTGACTCACCTCCACCAGACACGGACCACAGGGCCAACATGTGAGGATGAGCCG
ATCGTTACCTCAAGGGAGTGAAATTGGGCCAACTCCGGCCCTCCTGGTCACTTTTGGCCACGATGGC
CGGGGTACATACCTTGACCCGCGGAGGGCCAAGCGTAGTCCCAAGCATACCCACAGCGCTCCAGGAAGA
AGAATAAGAACTGCCGTGCCTTTCGCTCTACGTGGACTTCAGTGACGTGGGCTGGAATGATTGGATCGT
GGCCCCACCGGGTACCAGGCCTTCTACTGCCACGGGACTGTCCCTTTCCACTGGCGGACCACCTCAAC
TCAACCAATCATGCCATTGTGCAGACCCTGGTCAACTCCGTTAATTCTAGCATCCCTAAGGCCTGCTGTG
TCCCCACCGAAGTGAAGCCATTTCCATGTTGTATCTGGACGAGTACGACAAGGTGGTGTGAAAAATTA
TCAGGAGATGGTGGAGGGGTGCGGATGCCGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

ORF Size:	1224 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_012827.2 , NP_036959.2
RefSeq Size:	1559 bp
RefSeq ORF:	1227 bp
Locus ID:	25296
Cytogenetics:	15p14
MW:	46.5 kDa
Gene Summary:	plays a role in induction of cell proliferation [RGD, Feb 2006]