

Product datasheet for **RR212618**

Elavl1 (NM_001108848) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Elavl1 (NM_001108848) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Elavl1
Synonyms:	Hua; HuR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR212618 representing NM_001108848 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTAATGGTTATGAAGACCACATGGCGGAAGACTGCAGGGATGACATTGGGAGAACGAATTTAATTG
TCAACTACCTCCCTCAGAACATGACCCAGGAGAACTACGGAGTCTGTTCAGCAGCATTGGCGAGGTTGA
ATCTGCAAACCTATTTCGGGATAAAGTTGCAGGACACAGCTTGGGCTACGGCTTTGTGAAGTATGTTACT
GCAAAAGATGCAGAGAGCAATCAGCACACTGAACGGCTTGAGGCTTCAGTCCAAAACCATTAAGGTAT
CATATGCTCGCCCAAGCTCAGAGGTTATCAAAGATGCCAACTGTACATCAGTGGGCTTCCAAGGACCAT
GACTCAGAAGGATGTGGAAGACATGTTTTCTCGGTTTGGGCGAATCATCAACTCCAGGGTCTTGTGGAT
CAGACCACAGGTTTGTCCAGAGGGGTTGCCTTTATCCGGTTTGACAAACGGTCTGAAGCAGAAGAGGCAA
TTACCAGTTTCAATGGTCATAAACCCCCGGTTCCTCCGAGCCATCACAGTGAAGTTTGCAGCCATCC
CAACCAGAACAAAAACATGGCTCTCCTCTCGCAGCTGTACCACTCGCCTGCTAGGAGGTTTGGAGGCCCT
GTGCATCACCAGGCACAGAGATTCAGGTTCTCCCTATGGGTGTAGATCATAGAGTGGGATTTCTGGTG
TCAATGTCCCGGCAATGCTTCTCGGGCTGGTGCATCTTCATCTACAACCTTGGGCAAGACGCCGATGA
GGGATCTCTGGCAGATGTTTGGCCCCTTTGGTGCAGTTACCAATGTGAAAGTATTGTTGATTTCAAC
ACCAACAAGTGCAAAGGTTTGGTTTTGTGACCATGACAAACTATGAAGAAGCTGCAATGGCCATAGCAA
GTCTGAACGGCTACCGCCTGGGGGACAAAATTTACAGGTTTCTTCAAACCAACAAGTCCACAAAAGT
CTGAACGGCTACCGCCTGGGGGACAAAATTTACAGGTTTCTTCAAACCAACAAGTCCACAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR212618 representing NM_001108848
 Red=Cloning site Green=Tags(s)


MSNGYEDHMAEDCRDDIGRTNLI VNYLPQNMTQEELRSLFSSIGEVESAKLIRDKVAGHSLGYGFVNYVT
 AKDAERAISTLNLRLQSKTIKVSYPSEVIKIDANL YISGLPRTMTQKDVEDMFSRFGRINSRVLVD
 QTTGLSRGVAFIRFDKRSEAEAAITSFNGHKPPGSSEPI TVKFAANPNQKNMALLSQLYHSPARRFGGP
 VHHQAQRFRRFSPMGVDHMSGISGVNVPGNASSGWCIFIYNLQDDEGILWQMFPGPF GAVTNVKVIRDFN
 TNKCKGFGFVTMTNYEEAAMAIASLNGYRLGDKILQVSFKTNKSHKV*TATAWGTFYRFPSKPTSPT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1371_e11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_001108848

ORF Size: 978 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq Size:	2439 bp
RefSeq ORF:	981 bp
Locus ID:	363854
UniProt ID:	B5DF91
Cytogenetics:	12p12
MW:	36.2 kDa
Gene Summary:	RNA binding protein; involved in posttranscriptional regulation of VEGF expression [RGD, Feb 2006]

Product images:


Circular map for RR212618