

## Product datasheet for **RR212620L3V**

### Sh3gl3 (NM\_031238) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sh3gl3 (NM_031238) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sh3gl3
Synonyms:	Sh3d2c1; SH3P13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031238
ORF Size:	936 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212620).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031238.2</a> , <a href="#">NP_112517.1</a>
RefSeq Size:	1457 bp
RefSeq ORF:	939 bp
Locus ID:	81921
Cytogenetics:	1q31
Gene Summary:	interacts with inositol 5-phosphatase synaptojanin; may have a role in synaptic vesicle recycling [RGD, Feb 2006]



[View online »](#)