

Product datasheet for **RR212640L4V**

Ap2m1 (NM_053837) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ap2m1 (NM_053837) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ap2m1
Synonyms:	Ap50; Clapm1; mu2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053837
ORF Size:	1305 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212640).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053837.1
RefSeq Size:	1816 bp
RefSeq ORF:	1308 bp
Locus ID:	116563
UniProt ID:	P84092
Cytogenetics:	11q23
Gene Summary:	encodes the mu subunit of the heterotetrameric coat assembly protein complex 2 (AP2); associated with clathrin-coated vesicles and may be important for intracellular trafficking [RGD, Feb 2006]



[View online »](#)