

Product datasheet for **RR212666L3V**

Rfc2 (NM_053786) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rfc2 (NM_053786) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rfc2
Synonyms:	MGC95315
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053786
ORF Size:	1047 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212666).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053786.1 , NP_446238.1
RefSeq Size:	1556 bp
RefSeq ORF:	1050 bp
Locus ID:	116468
UniProt ID:	Q641W4
Cytogenetics:	12q12
Gene Summary:	40 kDa subunit of replication factor C; displays decreased expression in postnatal ventricular myocardium; may play a role in neonatal myocardial development and cell proliferation [RGD, Feb 2006]



[View online »](#)