

Product datasheet for **RR212685**

Clic3 (NM_001013080) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Clic3 (NM_001013080) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Clic3
Synonyms: MGC109225
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR212685 representing NM_001013080
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCTGAAACCACCAAGCTCCAGCTCTTTGTGAAGGCCAGTGAGGATGGCGAGAGCGTGGGACTGCC
CTTCCTGTGTCAGCGGCTTTCATGGTCTGCTCCTCAAGGGCGTGCCCTTCACACTACCACAGTGGACAC
ACGAAGGGCATTGGATGTGCTGAAGGACTTCGCACCGGGATCACAATTGCCATCCTACTCTATGATGGG
GACGTCAAGACAGACACTGCAGATTGAGGAGTTTCTGGAAGAGACATTGGGGCCACCCGACTTTCCTG
GCCTGGCACCTAGGTACAGGGAGTCCAACACAGCAGGCAATGACATCTTCCACAAGTTCTCTGCCTTCAT
CAAGAACCAGTGCCTACACAGGACGATGCCCTGTATCAGCAGCTGCTCCGGGCCCTGACCAGGCTGGAC
AGGTACCTGGGTACACACTGGACCATGAACTGGCACAGGAGCCACACCTTTCTGAATCCCGCCGCCGCT
TTCTCGATGGTGACCAGCTCACGTTGGCTGACTGCAGCCTGTACCCAAGCTGCACATTGTGGATACCGT
GTGTGCGCACTCCGCCAGAGGCCATCCCTGCGGAAGTGTCTGCGTTCCGCCGCTACTTGGACAGTGCA
TTGCAGGAGAAAAGTTCAAGTACACGTGTCCACACAGCGCAGAGATCCTGGCAGCCTACCAGCCTGCTG
TGCACCCCGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR212685 representing NM_001013080
 Red=Cloning site Green=Tags(s)

MAETTKLQLFVKASEDGESVGHCPSCQRLFMVLLKGVPTLTVDTRRALDVLKDFAPGSQLPILLYDG
 DVKTDTLQIEEFLEETLGPPDFPLAPRYRESNTAGNDIFHKFSAFIKNPVPTQDDALYQLLRALTRL
 RYLGTPLDHELAQEPHLSERRRFLDGDQLTLADCSLLPKLHIVDTVCAHFRQRP
 IPAE LSCVRRYLD SA LQEKEFKYTCPHSAEILAA YQPAVHPR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001013080
ORF Size: 711 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001013080.2 , NP_001013098.2
RefSeq Size:	792 bp
RefSeq ORF:	714 bp
Locus ID:	296566
Cytogenetics:	3p13
MW:	26.8 kDa