

Product datasheet for RR212686L4

Erg (NM_133397) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Erg (NM_133397) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Erg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212686).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

RefSeq ORF: 1368 bp
Locus ID: 170909
Cytogenetics: 11q11
Gene Summary: human homolog may act as a transcriptional activator [RGD, Feb 2006]