

Product datasheet for **RR212689**

Noc2l (NM_001033897) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Noc2l (NM_001033897) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Noc2l
Synonyms:	RGD1309387
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:


>RR212689 representing NM_001033897
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGATTTCTGTACTCTTTCCAGCTTTGTTCTTTGAATCTGGTACAGGAAACAAGTGAACCGGAGGAAG
ATGGAGGAGGGGACAGTGACACATTGCCAGAGGGCTGCAGAGGAAGAAGAATGAGCCTGTACCTGTGAC
CCTTGCCATGGTGGAGAGGTGGAAGCAGGGCGCCATGCACCACCTTACTCCCAAGCTGTTTCATGAAGTT
GTACAGGCATTCCGGGCAGCTGTAGCCACCACCAAGGAGACCAGGAAGCTCCTGAGACTTGCAGGTTCC
AGGTACAGACAGTGTGTGTTCAATGCTCTGGTTACCTTCTGCATTGAGACTCTGTGGTTGCCCTCA
GAAGCTACTGTTGCAAAGACAGTAAAGGATAGCAGCCGGATGCTGCTGCCCTCTAGTAGCCCACTGTGG
GGGAAGCTTCGTGTGGATGTGAAGACATATCTGAGTGCAGTTGTACAGCTGGTGCCTCTCTGGCTGAAG
CCACGGTGTCTGCCGCTGTCTGCAGCATATCAGCATTGTTGGTTCCTTACTACCTGACTTCCCTAAGCA
GTGCCGAATGCTGGTCAAGAGGATGATAGTTCTGTGGAGCACGGGTGAAGAATCTCTGCGGGTCTGGCC
TTCTGTGACTCATCAAAGTCTGCCGGCACAAGAAGGAATCCTTTCTTGGTCCCATCCTGAAGCAAATGT
ACATCATGTATGTGAGAACTGCAAGTTCACCTCTCCCAGTACCCTCCCCCTCATAAGCTTCATGCAGCG
GACACTGACGGAAGTCTTGCCCTGGACCCAGCGTCTCCTATCAGCACGCCCTTCTCTACATCCGACAG
CTTGACAGTCCACCTGCGTAATGCTATGACCACTGGCAAGAAGGAGACATAACAGTCTGTGTACAAGTGGC
AATACGTGCACTGCCTGTACCTGTGGTGTGAGTCTGAGTGCCTTGGTTCAGCGAGGTCCTGCAACC
CTTACTCTACCTCTCTCACAGATCATTTGGCTGTATCAAGTGTGTTGCCACTGCCCGCTTTTATCCG
TTACGCATGCATTGTGTACGTGCCCTGACACTGCTGTCCAGACCGTCCGTACCTTATCCCTGTTCTGC
CCTTCATTCTCGAGATTTTCCAACAGTGGACTTCAATAAGCGGCCAGGCCGATGAGCTCCAAGCCCAT
CAACTTCTCTGTGATCTTGAAGCTGTCCAGCACCAACCTGCAGGAGAAGGCATACCGGGATGGTCTGCTG
GAACAGCTGTATGACCTTATTATGGAATACCTGCACGGCCAGGCCACAGCATCGCTTCCCGGAGTTGG
TGCTGCCTACTGTTCTACAGCTGAAATCTTTTCTCCGGGAGTGCAAAGTGGCTAACTACTGCCGGCAGGT
GCGTCAGCTGCTGGAGAAAGTGAAGAGAACGCAGAACACATCCGAAGCCTTCGACAGAGAGTGACCTTC
AGTGTGTCTGACCAGCTGGCAGTGGATGCTTGGGAGAAGAAGTCCGTGAAGAGGGGACCCCTCTACCA
GATACTACAGCCACTGGAAGAAGCTGAGGGACCGTGAGATCCAGCTGGAGATCAGTGGCAAAGAGCGGT
AGAAGACCTGAAGTCCCAGAGATCAAAGGCGGAAAGTGAAGACAGGAAGGATGAAGACAGGAAAGAA
TTAAAAGACCTGTTTGGATTGGACAGTTCTGAGGGAGAGGACAGCACAGACTTCTTTGAGAGAGGAGTTC
CTGGGCTCCTGGAATCTCACAGAGGAATGAAAGAAGGAGGTGAGGAGGAAGAACAATAAGGAAGATAGTGA
CAGCAGTTCAGAGGATGGAGACACAGACACAGGGTGGATCTGAGTGAAGTGTGGCAGCTGGCTCAGGGA
CCACAGGATGAGCTAGAGGATCTTCAGCTCTCAGAAGAGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:



ACCN: NM_001033897

ORF Size: 1932 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033897.1](#), [NP_001029069.1](#)

RefSeq Size: 2485 bp
RefSeq ORF: 1935 bp
Locus ID: 313777
Cytogenetics: 5q36
MW: 73.8 kDa