

Product datasheet for **RR212692L3V**

Akap8 (NM_053855) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Akap8 (NM_053855) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Akap8
Synonyms:	Akap95
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053855
ORF Size:	2061 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212692).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053855.1 , NP_446307.1
RefSeq Size:	3438 bp
RefSeq ORF:	2064 bp
Locus ID:	116633
UniProt ID:	Q63014
Cytogenetics:	7q11
Gene Summary:	scaffolding protein that is involved in targeting type II cAMP-dependent protein kinase for cAMP-responsive nuclear events [RGD, Feb 2006]



[View online »](#)