

Product datasheet for **RR212743L4V**

Bet1 (NM_019251) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bet1 (NM_019251) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Bet1
Synonyms:	MGC72488
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019251
ORF Size:	354 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212743).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019251.2 , NP_062124.1
RefSeq Size:	1455 bp
RefSeq ORF:	357 bp
Locus ID:	29631
UniProt ID:	Q62896
Cytogenetics:	4q13
Gene Summary:	protein involved with vesicular membrane traffic, especially with the golgi; may be associated with the SNARE complex [RGD, Feb 2006]



[View online »](#)