

## Product datasheet for **RR212758**

### Blvra (NM\_053850) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Blvra (NM\_053850) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Blvra  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR212758 representing NM\_053850  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGATGCCGAGCCAAGAGGAAATTTGGAGTGGTAGTGGTTGGTGTGGCAGAGCTGGCTCGGTGAGGC  
TGAGGGACTTGAAGGATCCACGCTCTGCAGCATTCCTGAACCTGATTGGATTTGTGTCCAGACGAGAGCT  
TGGGAGCCTTGATGAAGTACGGCAGATTTCTTTGGAAGATGCTCTCCGAAGCCAAGAGATTGATGTCGCC  
TATATTTGCAGTGAGAGTTCCAGCCATGAAGACTATATACGGCAGTTTCTGCAGGCTGGCAAGCATGTCC  
TCGTGGAATACCCCATGACACTGTCATTTGCGGGCCAGGAGCTGTGGGAGCTGGCCGCACAGAAAGG  
GAGAGTCTGCATGAGGAGCACGTGGAACCTTTGATGGAGGAATTCGAATTCCTGAGAAGAGAAGTGTG  
GGGAAAGAGCTACTGAAAGGGTCTCTTCGCTTACAGCTAGCCACTGGAAGAAGAGAGATTTGGCTTCC  
CTGCGTTCAGCGGCATTTCTCGCCTGACCTGGCTGGTCTCCCTCTTCGGGGAGCTTTCTTATTCTGC  
CACCTTGAAGAGCGAAAAGAGGATCAGTATATGAAAATGACCGTGCAGCTGGAGACCCAGAACAAGGT  
CTGCTGCATGGATTGAAGAGAAAGGGCCTGGCTTAAAAAGAAACAGATATGAAACTTCCAGTTCACTT  
CTGGGTCCCTGGAGGAAGTGCCAAGTGTAGGGGTCAATAAGAACATTTTCTGAAAGATCAGGATATATT  
TGTTCAGAAGCTCTTAGACCAGGTCTCTGCAGAGGACCTGGCTGCTGAGAAGAAGCGCATCATGCATTGC  
CTGGGGCTGGCCAGCGACATCCAGAAGCTTTGCCACCAGAAGAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR212758 representing NM\_053850  
 Red=Cloning site Green=Tags(s)

MDAEPKRKFGVVVVGVRAGSVRLRDLKDPRSA AFLNLIGFVSRR ELGSLDEV RQISLEDALRSQEIDVA  
 YICSESSSHEDYIRQFLQAGKHV LVEY P M T L S F A A A Q E L W E L A A Q K G R V L H E E H V E L L M E E F F L R R E V L  
 GKELLKGS LRFTASPLEEERFGFP AFSGI SRLTWLVSLFGELSLISATLEERKEDQYMKMTVQLETQNKG  
 LLSWIEEKGPGLKRNRVNFQFTSGSLEEVPSVGVNKNIFLKDQDIFVQKLLDQVSAEDLAAEKKRIMHC  
 LGLASDIQKLCHQKK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_053850

<b>ORF Size:</b>	885 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_053850.1</a> , <a href="#">NP_446302.1</a>
<b>RefSeq Size:</b>	1082 bp
<b>RefSeq ORF:</b>	888 bp
<b>Locus ID:</b>	116599
<b>UniProt ID:</b>	<a href="#">P46844</a>
<b>Cytogenetics:</b>	3q36
<b>MW:</b>	33.6 kDa
<b>Gene Summary:</b>	catalyzes the conversion of biliverdin to bilirubin in heme degradation [RGD, Feb 2006]