

Product datasheet for RR212764L4

Cnot6l (NM_001108355) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cnot6l (NM_001108355) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Cnot6l
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212764).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001108355
ORF Size:	1650 bp



[View online »](#)

