

Product datasheet for **RR212786**

Tnfaip6 (NM_053382) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tnfaip6 (NM_053382) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tnfaip6
Synonyms:	Tnfip6
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR212786 representing NM_053382 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGTCTTCCTCCTTTGCTTATGCGTCTTGCTGTGGGAAGAGGCTCACGGATGGGGATTCAAGAATGGAA
TCTTTTCATAATCCATATGGCTTGAACAAGCAGCAGGCGTATACCATAGAGAAGCCCGGGCTGCCAGATA
CAAGCTAACCTATGCAGAAGCCAAGGCGGTATGTGAATATGAAGGTGGCCGTCTTCAACCTACAAGCAG
CTAGAAGCAGCCAGAAAGATTGGATTCCATGTCTGTGCTGCTGGGTGGATGGCCAAGGGTAGAGTTGGGT
ACCCCATTTGTGAAACCTGGGCCAACTGTGGATTTGGCAAAACCGGCATCATTGATTATGGAATCCGGCT
CAACAGGAGCGAGCGGTGGGACGCCTATTGCTACAACCCGAACGCAAAGGAGTGTGGTGGTGTCTTTACA
GACCCGAAGCGAATCTTTAAATCCCGAGGCTTCCCAAATGAGTATGACGACAACCAGGTTTGCTACTGGC
ACATCCGGCTCAAGTACGGTCAGCGAATCACCTGAGCTTTTGGACTTTGACCTTGAACATGATCCAGG
CTGTTTGGCTGACTATGTAGAAATATGACAGTTATGACGATGTCCACGGCTTTGTAGGAAGATACTGC
GGTATGAACTTCCAGAAGACATCATCAGCACAGGAAATGTCATGACCTTGAAGTTTCTGAGTATGCGT
CCGTCACAGCCGGAGGCTTCCAGATTAATACGTCACAGTGGATCCTGCACCTAAATCCAATCAAGGCCAA
AAATACAACCTATGACTGGAATAAAAAGTTCTTAGCTGGACGGTTTAGCCATCTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR212786 representing NM_053382
Red=Cloning site Green=Tags(s)

MVFLLCVLLWEEAHGWGFKNGIFHNSIWLEQAAGVYHREARAGRYKLTAEAKAVCEYEGRLATYKQ
 LEAARKIGFHVCAAGWMAKGRVGYPIVKPGPNCGFGKTGIIDYGIRLNRSERWDAYCYNPNAKECGGVFT
 DPKRIFKSPGFPNEYDDNQVCYWHIRLKYGQRIHLSFLDFDLEHDPGCLADYVEIYDSYDDVHGFVGRYC
 GDELPEDIISTGNVMTLKFLLSDASVTAGGFQIKYVTVDPAPKSNQGKNTTMTGNKKFLAGRFSHL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_053382
ORF Size: 825 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_053382.1 , NP_445834.1
RefSeq Size:	1401 bp
RefSeq ORF:	828 bp
Locus ID:	84397
Cytogenetics:	3q12
MW:	31 kDa
Gene Summary:	mouse homolog is a catalyst in the formation of the cumulus extracellular matrix and indispensable for female fertility [RGD, Feb 2006]