

Product datasheet for **RR212819**

Dlat (NM_031025) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dlat (NM_031025) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dlat
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR212819 representing NM_031025
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTGGCGCTGTGTGCGCGGAGGGTGCAGAGTGTCTGCCAGGGCCGGTTTCAGGGCTCGGTGGGCAA
 CTTTGAAGGGCCCGCGGACCGCCCGCCGAGTTCGTTGCGGCAGTGGGATACCCAGCTATGGGGTCCG
 CTGTTGTGCGGCTGGAGCTACGGCTCTGCCACGTCACCGTAACCGCATACTGCAACAACCTCCTGGG
 TCGCCAGCCGCGCTCTATAGCCTCCCTCCGCACCAGAAGGTTCCACTACCTTCACTTTCCCTACAA
 TGCAGGCAGGCACCATAGCCCGTGGGAAAAAAGAAGGGAAAAAATCAGTGAGGGTGTCTAATTGC
 AGAGGTTGAAACTGATAAAGCCACTGTTGGATTTGAGAGCTTGGAGGAGTGTATATGGCAAAGATCCTT
 GTTCTGAGGGGACCAGAGATGTTCCAGTAGGATCTATCATTTGTATCACAGTTGAAAAGCCTCAGGATA
 TTGAGGCCTTAAAAATTATACATTGGATTGAGCAACAGCTGCCACACAGGCAGCCCGAGCCCGAGCTGC
 AGCTCCAGCTGCTGCCCGAGCTGCACCTTCTGCAAGTGTCCAGGTAGCTCCTATCCCGTTCACATGCAG
 ATTGTTCTTCTGCCCTCTCCCAACCATGACCATGGGCACCGTCCAGAGGTGGGAAAAAGAAAGTGGGAG
 AGAAGCTGAGTGAAGGAGACTTGTCTGGCAGAGATAGAGACCGACAAGGCCACCATAGGCTTTGAAGTACA
 AGAAGAAGGTTATCTGGCAAAAATCCTGGTCCCTGAAGGCACAAGGGATGTTCTCTGGGAACCCCGCTA
 TGTATCATAGTAGAAAAACAGGAAGATATAGCAGCATTTCGAGACTACAGGCCAACAGAAGTGACCAGCT
 TAAAGCCACAGGCACCACCACCTGTCCCACCCAGTGGCAGCTGTTCTCCCATCCCCAGCCTTTAGC
 ACCTACCCCTCAGCCGCTCCTGCTGGACCAAGGGAAAGGGTGTTCGTTAGCCCTCTTGCAGAAAGATTG
 GCAGCAGAGAAAGGATTGACCTCACCAAGTTAAAGGGACGGGACCAGAAGGCAGAATCATCAAGAAGG
 ACATTGACTCTTTTGTGCTACTAAGGCTGTCTGCGCTGCAGCAGCTGCTCCCCGGTCCAAGAGT
 GGCACCAACTCCTGCAGGTGTCTTCATAGACATCCCATCAGCAACATTGTCGAGTGATTGCGCAGAGG
 CTCATGCAGTGAAGCAGACTATACCTCATTATTACCTTTCTGTTGATGTAATATGGGAGAGGTGCTGT
 TGGTGGGAAGGAACTTAATAAGATGCTTGAAGGTAAAGGAAAAATCTCCGTCAATGACTTCATCATAAA
 AGCTTCAGCTTTGGCCTGCCTGAAAGTTCCTGAAGCAAACCTCATCTTGGATGGACACAGTTATACGACAA
 AATCATGTGGTTGATGTCAGCGTTGCTGTCAGTACCCCTGCAGGACTTATCACCCCTATTGTGTTAATG
 CACACATAAAAGGACTGGAAACCATTGCTAGTGATGTTGTTTCTTAGCCTCAAAGCAAGAGAGGGTAA
 ACTTCAGCCTCACGAGTTCAGGGTGGGACATTTACAATCTCAAACCTTAGGGATGTTCCGAATTAAGAAT
 TTCTCTGCGATTATTAACCCACCTCAGGCATGATTTTGGCAATTGGTGTCTCCGAGGATAAACTGATCC
 CTGCAGATAATGAGAAAGGCTTTGACGTGGCTAGTGTGATGTCTGTACACTCAGCTGTGATCATCGAGT
 TGTGGATGGAGCAGTTGGAGCCAGTGGCTTGCTGAATTCAAAAAGTACCTTGAAAAGCCTGTAACATATG
 TTATTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR212819 representing NM_031025
 Red=Cloning site Green=Tags(s)


MWRVCARRVQSAVPRAGFRARWATLKGPRTPGAAVRCGSGIPSYGVRSLCGWSYGSATVPRNRILQQLLGS
 SPSRRSYSLPPHQKVLPSLSPTMQAGTIARWEKKEGEKISEGDLIAEVETDKATVGFESLEECYMAKIL
 VPEGTRDVPVGSIIICITVEKPQDIEAFKNYTLDSATAATQAAPAAAAAAPAAPSASAPGSSYPVHMQ
 IVLPALSPTMTMGTVQRWEKKVGEKLESEGDLAEIETDKATIGFEVQEEGYLAKILVPEGTRDVPVPLG
 CIIVEKQEDIAAFADYRPTVETSLKPQAPPPVPPVAAVPPVPPVAAVPPVPPVAAVPPVPPVAAVPPV
 AAEKIDLTQVKGTGPEGRIKKDIDSFVPTKAAAPAAAAAAPPGRVAPTPAGVFDIPIISNIRRVIAQR
 LMQSKQTIPIHYLLSVDVNMGEVLLVRKELNKMLEGKGISVNDFIKASALACLKVEANSSWMDTVIRQ
 NHVVDVSVAVSTPAGLITPIVFNHAIKLETIASDVVSLASKAREGKLQPHFQGGTFTISNLGMFGIKN
 FSAIINPPQACILAIIGASEDKLIPADNEKGFVDVSVMSVTLSCDHRVVDGAVGAQWLAEFKKYLEKPVMT
 LL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:


SgfI-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_031025

ORF Size: 1896 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_031025.1](#), [NP_112287.1](#)
RefSeq Size: 2628 bp
RefSeq ORF: 1899 bp
Locus ID: 81654
UniProt ID: [P08461](#)
Cytogenetics: 8q23
MW: 67.2 kDa
Gene Summary: putative dihydrolipoamide acetyltransferase [RGD, Feb 2006]