

Product datasheet for **RR212874**

Tmem268 (NM_001106661) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmem268 (NM_001106661) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tmem268
Synonyms: RGD1304595
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR212874 representing NM_001106661
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTTGTAACCTCCCTCAGATCCTGGTGGGGCAGCTGGTCTCTGCCACCTCCACCCTGGGCTGCA
GCACCCTGCCTCAGGGGAACCTCCTGGCTGGGGGAAGAGCTACACAATGGCCAGTCTCTCAGTCTCT
CCGGATCGACAACACCTGTGCACCCATCTCCTTTGACCTGAGAGTTGCAGAAGAAGAGCTGCAGGCCTGG
GGCATTACAGTCCCTGCTGAGCAGTATAGGAGCTTGGCTGAGAGTGCCTTTTGGAGCCTCAAGTGAGGA
GGTACATCATCTACAACCTCCAGGCCATGAGGCTGGCCTTTGCTGTGGTGTCTATGTGTTGGTGTGGGC
CAACATCTACTCCACCAGCCAAATGTTTGCCCTAGGCAACCAGTGGGCAGGTGTGCTGCTTGCACCCCTG
GCCGCTGTGAGCCTGACCCTGACTCTTGTGCTGGTCTTTGAGAGGCAGCAGAGGAAGGCTAACACTAACA
CAGACCTGAGGCTGGTGGCTGCCAATGGAGCCCTTCTGAGACACCGAGTGTGCTGGGGGTGACAGATAC
CGTGGAGGGATGCCAGAGTGTGATACAGCTCTGGTTTGTCTACTTTGACCTGGAGAAGTGTGTGCAGTTC
TTGTCCGATCATGTCCAAGAAAATGAAGAGGAGTCAGGAGTCCCTTGTGAGAAGCAGATTAAGCCAGCTAT
GTGTTGTCATGGAGACTGGAGTGGAGTGGCCCGTGGTGGAAAGGCCTGAGGACTTGGAAAGATCTCCCCTCT
GCCCAGAACTCCTGGTCCCTCAGGAGAGACCACTCACACAGACTGAACTCTATCAGCTTGTCCAGAGGCT
GAGCCAGAGGAGATGGCCCGGCAACTGCTGGCAGTGTGGTGGTACTACTACCCCGGCTCCTGGTACCT
CTCAGCTCCCTCAGCCAATGGGAACACGACACATGGACTCTGCAAGGATCCCATGCCCTGCCAGCTCAT
AGAAGTACACATTTTGGGCACAGGATGTTGCCATTCTGGCCAGG

ACGGTACGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001106661.1 , NP_001100131.1 |
| RefSeq Size: | 3795 bp |
| RefSeq ORF: | 1029 bp |
| Locus ID: | 298104 |
| Cytogenetics: | 5q24 |
| MW: | 37.9 kDa |