

Product datasheet for **RR212930L3V**

Nostrin (NM_001024260) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nostrin (NM_001024260) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nostrin
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001024260
ORF Size:	1506 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212930).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024260.1 , NP_001019431.1
RefSeq Size:	1720 bp
RefSeq ORF:	1509 bp
Locus ID:	311111
UniProt ID:	Q5I0D6
Cytogenetics:	3q21
Gene Summary:	Multivalent adapter protein which may decrease NOS3 activity by inducing its translocation away from the plasma membrane.[UniProtKB/Swiss-Prot Function]



[View online »](#)