

## Product datasheet for **RR212966**

### Clca4 (NM\_201419) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Clca4 (NM_201419) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Clca4
Synonyms:	Clca6; Prp3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RR212966 representing NM\_201419  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGATGGGGTTCTCCAGAGGGATTGTTTTCTCTCTGCTCTGTACCTGCTTCAGGGGTCAGACACCTCT  
 TGGTAAAGCTGAATGAGAATGGTTATGAAGACATTATTATTGCCATCGATCCTGCTGTGTCCGAAGATGT  
 AACAAATATTGACCAAATAAAGGACATGGTGACTAAAGCATCTGCATACCTTTTTGAAGTACAGAAAA  
 AGATTTTTTTTTCAAAACGTCTCTATATTAATTCCTGAGAATTGGACAAACAGTGACCAGTACAGGAGGC  
 CAAAGCAGGAAAGCTACAACATGCTGACATTAAGTTGCACCTCCTGCTCTCCAGGGCAGGGATGAGCC  
 CTACACCCGGCAGTTCACAAAGTGTGAAAGAAAGCTGAATACATCCACTTACCCCCGACTTTGTATTG  
 GGGAGAAAAACAAAGGAATATGGAGATTGAGAAAGACTGTTGGTCCATGAATGGGCCATCTTCGCTGGG  
 GAGTATTTGATGAGTACAACGAAGACCAGCCATTCTACAGTGCTTCATCAAAGAAAAATCGAAGCAACAAG  
 GTGTTCCACAGGTATCAAGGGTATGAATAAAGCCAGGTATGCCAAGGGGGCAGCTGTATAACCAGAAAT  
 TGCAGACGGAATTCACCACACAGTTATATGAAAAAGATTGTCAGTTCTCCCGGATAAAGTTCAAACAG  
 AAAAGTCATCTATAATGTTTCATGCAAAGTATCGATTCACTGACTGAATTCGTAAAAAGAAAAACATAA  
 TCGAGAGGCTCCAACCTACATAATCAGAAGTGCATTACAGAAGCACATGGGAAGTAAATCAGTAATTCT  
 GAGGATTTTTAAAAACAGCACACCCATGGAGATGCCACCCTCTCCACCCTTCTTCGCTGCTGAGGATCA  
 GTGAAAGAATTGTGTGCTAGTGCTTGTGTCTGGGAGCATGGGGTCTTATGATCGCCTCAATCGAAT  
 GAACCAAGCAGCAAAAATCTTCTACAGCAGATACTTGAAGTAGGTCTGGGCAGGGATGGTACACTTT  
 CATAGCTCGGCCACTGTTAAGAGTGAGCTGATCCAAAATAACAGTGATGTCGAGAGAAACCAGCTTTTGG  
 AAACCTTACCTACAAGTGTAGCGGAGGGACATCCATCTGCTCTGGAATCAGAAGTCAATTCAGGTATT  
 TAAAAATAAAGGATACCAAACCTGGTGGGAATGACATCTTGTGTTGAGCGATGGGGAGGACAGCACAGCC  
 AAAGACTGCCTTGATGAGGTGAAAGACAGCGGGCCGTAGTTCATTTTCATTGCCTTGGGAAAAGCATTG  
 ATCAATCTATATCAAACATGGCCAACGTAACAGGTGGAAGCAACTGTTTGTACAGATGAAGCCAGAA  
 CAATGGACTCATTGATGCTTTCGGGGCCCTGGCTTCAGAAAAATGCTGATGTCAGTGAAGTCCCTTCAG  
 CTTGAAAGCAAGGGTGTGTGTGACCAACAGTCGCTGGCTGAACGACTGTCGTTATTGACAGCACCG  
 TGGGAAAGGACACATACTTTCTGTACCTGGAGCCAACAAGCCCTGCCATTCTCAGGGATCCCAA  
 AGGAACACAAATAACAAATTTACAATCGACACGGCTTCAAAAATGGCCTATCTCACTATCCCAGGAATT  
 GCTGAGGTGCGGCTTTGACTTACAACCTTGAAGCCAAAGTAGACTCAGAGATACTAACAATTACAGTAA  
 CTTCTCGGCATCACGCTCTTGTGCCACCCATCATTGTAATGCAAAAAGTTAACACAGACACGAAACAG  
 TTTCCCAAGTCCAATGATTGTGTATGCAGAAGTTCTACAAGTTACACTCCCATCATTGGAGCACGCGTG  
 ACAGCTACCATAGAATCAAAATAGTGGAAAAACAGAAGAGCTAGTACTGCTGGACAATGGTGCAGGGGCTG  
 ATGCTTCAAGGATGACGGGGTCTACTCCAGGTACTTCACGGCTTATTAGATAATGGCAGATATAGCGT  
 AAAGGTCCGCGCAGACGGAGGAACAACTCTGCCCGGAGAAGTTCACGGCATCCGTCCAGCAGAGCGGCA  
 TACATACCAGGCTGGGTAGTGGATGGAGAAATCAAGGGAACCCACCCAGACCTGAAACTACTGAGGACA  
 CTCAGCCCGTCTGGAGAATTTACAGCAGGACGGCATCTGGAGGTGCCTTTGTGATGTCCAAGTCCCAGC  
 TCTTCCATTGCCTGATCAGTACCCACCAATCAAAATCACAGACCTCCAGGCCACACTTGATGGGAAAAAG  
 ATCAGTCTAACATGGACGGCACCAGGAGATGACTTTGATGTTGGAAGAGTTCAACAGTATATCATAAGAA  
 CAAGTAAAAATATCATTGACCTCAGAGACAATTTAATAAGTCTCTTCAGGTTAATACAACTCACCTGAT  
 ACCAAAAGAAGCCAACCTCAAGGAAACCTTTGCCTTTAAACCAGAAACTATCTCAGGAGAAAAACGCAACC  
 TACATATTCAATGCCATTGAAAGTGTGACAAAAAGCAATTTGCACTCCGGACTGTCCAACATCGCACAAG  
 TAGCACTGTTACCCCTCAGGCAGATCCTGTCCCTGAGGAAAGTCCCAGTCTGGAGTTAGTATCACTAC  
 GATTGTGCTGTGTGGTGGGCGCGTGTACTAGTTTGTATCATTGTAAGTACCAGTCTTTGTATCTTA  
 AAGAGAAAAAGGTCCTCCTCAGGGCCTGCAACAACATTT

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR212966 representing NM\_201419  
 Red=Cloning site Green=Tags(s)


MMGFSRGIVFLLLL YLLQGS DTS L V K L N E N G Y E D I I I A I D P A V S E D V T I I D Q I K D M V T K A S A Y L F E A T E K  
 R F F F K N V S I L I P E N W T N S D Q Y R R P K Q E S Y K H A D I K V A P P A L Q G R D E P Y T R Q F T K C G K K A E Y I H F T P D F V L  
 G R K Q K E Y G D S G R L L V H E W A H L R W G V F D E Y N E D Q P F Y S A S S K K I E A T R C S T G I K G M N K A Q V C Q G G S C I T R N  
 C R R N S T T Q L Y E K D C Q F F P D K V Q T E K S S I M F M Q S I D S V T E F C K K E N H N R E A P T L H N Q K C D Y R S T W E I S N S  
 E D F K N S T P M E M P P S P P F S L L R I S E R I V C L V L D V S G S M G S Y D R L N R M N Q A A K F F L Q Q I L E S R S W A G M V H F  
 H S S A T V K S E L I Q I N S D V E R N Q L L E T L P T S A S G G T S I C S G I R T A F Q V F K N K G Y Q T G N D I L L L S D G E D S T A  
 K D C L D E V K D S G A V V H F I A L G K A F D Q S I S N M A N V T G G K Q L F A T D E A Q N N G L I D A F G A L A S E N A D V T E K S L Q  
 L E S K G A V L T N S R W L N D T V V I D S T V G K D Y F L V T W S Q Q A P A I H L R D P K G T Q I T N F T I D T A S K M A Y L T I P G I  
 A E V G V W T Y N L E A K V D S E I L T I T V T S R A S R P S V P P I I V N A K V N T D T N S F P S P M I V Y A E V L Q G Y T P I I G A R V  
 T A T I E S N S G K T E E L V L L D N G A G A D A F K D D G V Y S R Y F T A Y S D N G R Y S V K V R A D G G T N S A R R S S R H P S S R A A  
 Y I P G W V D G E I Q G N P P R P E T T E D T Q P V L E N F S R T A S G G A F V M S K V P A L P L P D Q Y P P N Q I T D L Q A T L D G E K  
 I S L T W T A P G D D F D V G R V Q Q Y I I R T S E N I I D L R D N F N K S L Q V N T H L I P K E A N S K E T F A F K P E T I S G E N A T  
 Y I F I A I E S V D K S N L H S G L S N I A Q V A L F T P Q A D P V P E E S P Q S G V S I T T I V L S V V G A V V L V C I I V S T T V C I L  
 K R K R S S S G P A T T F


SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-RsrII

**Cloning Scheme:**



**Plasmid Map:**


**ACCN:** NM\_201419

**ORF Size:** 2769 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_201419.1](#), [NP\\_958822.1](#)

**RefSeq Size:** 3097 bp

**RefSeq ORF:** 2772 bp

**Locus ID:** 362053

**Cytogenetics:** 2q44

**MW:** 101.9 kDa