

Product datasheet for **RR212985**

Zbtb26 (NM_001107840) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zbtb26 (NM_001107840) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Zbtb26
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR212985 representing NM_001107840 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCCCTTCCCGGGTTGTCTGCCAAATGTCTGAAAGATCAGATCTCCTTCACTTCAAGTTT
AAAATTATGGAGATTCAATGTTACAAAAATGAACAAATTAAGAGAAGAGAATAAAATTTGTGATGTTAC
AGTTCTCATAGATGATGTTGAGGTGCAGGGGCATAAAATGTTTTGCTGCAGGTTCCCGCTTCTTAAGA
GACCAATTTTACTGAATGATCCCGAGAGGTGAAAATCTCCATATTACAGAGTTCTGAAGTGGGGAGAC
AATTACTTATCCTGTTATAGTGGTGTCTGGAATCCCTGAGATGGAAGTAAATTAAGTTGACTGCTG
TGCAAGTTTTCTCAGATGAGCCACATTGTAGAACGGTGCACACAGGCGCTGTGGAAGTTATAAAGCCA
AAACAACCAATGGATAGTAAAGAGGGATCTGAACCACAGAGTGTCTCCCGAGTCGAAAGAACAGCAGG
GAGATGCCAGAGGCTCTCCCAAACAGGACTTACCTTGTATTATCCATCTGAAGACAGTATGGATATGGA
GGACAGTGATATTCAGATTGTTAAGGTAGAATCTATTGGGGATGTATCAGAGGATAGAAGTAAAAAAGAT
CAGAACCAGTACATTTCTTGAACCCACTGCTTACATTCATCAGAGCCCCAGCACTCCCTGATAAACT
CAACTGTGAAAAACAGAGTTAATGAACTAGACCAGAGCCATCTCCACAATTACGCCCTCTTTACACAGG
CAATGATGACATCATCATGACCTCAAAGATGTCTTTGGGCCATAATTCGAGGTGTAGACAAAGGCTTA
CAGTGGCATCACCAATGCCCAAAGTGTACCAGGGTATTCCGTCACCTGGGAACTACGCCAACATTTAA
AAATGCACAACTCTTTATGTGTCTACTCTGCGGCAAGACTTTCAGTCAAGAAAGGCAATCTTCATCGACA
CATGCGTGTACATGCCGGAATTAACCTTTCCAGTGTAAAATCTGTGGGAAAACCTTTTCTCAGAAGTGT
TCCTTACAGGACCATCTAACCTTACAGTGGAGATAAGCCCCATAAATGTAAGTATTGTGACATGGTTT
TTGCACATAAGCCAGTGTGAGGAAACACCTTAAACAGTTGCATGGCAAAAACAGCTTTGATAATGCCAA
TGAGAGGAGTGTGCAAGACCTCACAGTGGATTTGACTCTTTGCATGTACAACAGTACAGACTCTAAT
TGTCAGCCACAACCTGATGCAACACAGGTCCTGGATGCAGGTAAGTACCCAGCTGTCTCAGCTTAA
GGAGTGACAGTACATGTGTAAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR212985 representing NM_001107840
Red=Cloning site Green=Tags(s)

MAAFPPGLSAKMSERSDLLHFKFENYGDSMLQKMNKLREENKFCDVTVLIDDDVEVQGHKIVFAAGSPFLR
 DQFLLNDSREVKISILQSSEVGRQLLLSCYSGVLEFPPEMELVNYLTAASFLQMSHIVERCTQALWKFIPK
 KQPMDSKEGSEPPQSASPQSKEQQGDARGSPKQDLPCIHPSSEDSMDMEDSDIQIVKVESIGDVSSEDRSKKD
 QNQYISSEPTALHSSEPPQHSLINSTVENRVNELDQSHLHNYALSYTGNDIDIIMTSKDVFGPNIRGVDKGL
 QWHHQCPKCTRVFRHLENYANHLKMHKLFMCLLCGKTFQKGNLHRHMRVHAGIKPFQCKICGKTFQK
 SLQDHLNLHSGDKPHKCNVCDMVFAHKPVLRKHLKQLHGKNSFDNANERSVQDLTVDFDSFACTTVTDSN
 CQPQPDATQVLDAGKLTQAVLSLRSDSTCVN


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001107840
ORF Size:	1353 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107840.1 , NP_001101310.1
RefSeq Size:	2219 bp
RefSeq ORF:	1356 bp
Locus ID:	311910
Cytogenetics:	3q11
MW:	51 kDa