

## Product datasheet for **RR213002L3V**

### Chn1 (NM\_032083) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Chn1 (NM_032083) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Chn1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032083
ORF Size:	1002 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213002).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032083.1</a> , <a href="#">NP_114472.1</a>
RefSeq Size:	2309 bp
RefSeq ORF:	1005 bp
Locus ID:	84030
UniProt ID:	<a href="#">P30337</a>
Cytogenetics:	3q23
Gene Summary:	expressed only in neurons; may mediate neuronal signal transductions; human homolog functions as GTPase-activating protein (GAP) for p21rac and phorbol ester receptor [RGD, Feb 2006]



[View online »](#)