

Product datasheet for **RR213028L4V**

Ptpn20 (NM_001080149) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ptpn20 (NM_001080149) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ptpn20
Synonyms:	Ptpn20b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001080149
ORF Size:	1263 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213028).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001080149.1 , NP_001073618.1
RefSeq Size:	3040 bp
RefSeq ORF:	1266 bp
Locus ID:	306281
UniProt ID:	A1L1L3
Cytogenetics:	16p16
Gene Summary:	Tyrosine-protein phosphatase targeted to sites of actin polymerization in response of varied extracellular stimuli. Has tyrosine phosphatase activity towards various tyrosyl phosphorylated substrates (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)