

Product datasheet for **RR213084L4V**

Adad1 (NM_001006977) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adad1 (NM_001006977) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adad1
Synonyms:	Tenr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001006977
ORF Size:	1857 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213084).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001006977.2 , NP_001006978.2
RefSeq Size:	1944 bp
RefSeq ORF:	1860 bp
Locus ID:	294979
UniProt ID:	Q3KR54
Cytogenetics:	2q25
Gene Summary:	Plays a role in spermatogenesis. Binds to RNA but not to DNA (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)