

Product datasheet for **RR213090L4V**

Pecam1 (NM_031591) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pecam1 (NM_031591) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pecam1
Synonyms:	CD31; Pecam
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031591
ORF Size:	2034 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213090).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031591.1 , NP_113779.1
RefSeq Size:	3244 bp
RefSeq ORF:	2037 bp
Locus ID:	29583
UniProt ID:	Q3SWT0
Cytogenetics:	10q32.1
Gene Summary:	cell adhesion molecule, also intracellular signaling molecule; may mediate leukocyte adhesion and migration, angiogenesis, and thrombosis [RGD, Feb 2006]



[View online »](#)