

## Product datasheet for **RR213092L4V**

### Usp30 (NM\_001107153) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Usp30 (NM_001107153) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Usp30
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001107153
ORF Size:	1551 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213092).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001107153.1</a> , <a href="#">NP_001100623.1</a>
RefSeq Size:	3362 bp
RefSeq ORF:	1554 bp
Locus ID:	304579
UniProt ID:	<a href="#">D3ZPG5</a>
Cytogenetics:	12q16



[View online »](#)

**Gene Summary:**

Deubiquitinating enzyme tethered to the mitochondrial outer membrane that acts as a key inhibitor of mitophagy by counteracting the action of parkin (PRKN): hydrolyzes ubiquitin attached by parkin on target proteins, such as RHOT1/MIRO1 and TOMM20, thereby blocking parkin's ability to drive mitophagy (PubMed:24896179). Preferentially cleaves 'Lys-6'- and 'Lys-11'-linked polyubiquitin chains, 2 types of linkage that participate to mitophagic signaling. Does not cleave efficiently polyubiquitin phosphorylated at 'Ser-65' (By similarity). Acts as negative regulator of mitochondrial fusion by mediating deubiquitination of MFN1 and MFN2 (By similarity).[UniProtKB/Swiss-Prot Function]